



PROTEIN
BIOTECHNOLOGIES

HUMAN UTERINE CERVIX CANCER TISSUE MICROARRAY

Catalog Number: TMA-005

Each core represents one specimen selected and pathologically confirmed. Normal matched or unmatched uterine cervix tissue is provided for comparison to the cancer tissue.

Cases: 68 **Cores:** 80 **Diameter:** 1.5 mm **Thickness:** 5 µm **Standard IHC:** Cytokeratin confirmed

Layout: 8 x 10

	1	2	3	4	5	6	7	8	9	10
A	1	2	3	4	5	6	7	8	9	10
B	11	12	13	14	15	16	17	18	19	20
C	21	22	23	24	25	26	27	28	29	30
D	31	32	33	34	35	36	37	38	39	40
E	41	42	43	44	45	46	47	48	49	50
F	51	52	53	54	55	56	57	58	59	60
G	61	62	63	64	65	66	67	68	69	70
H	71	72	73	74	75	76	77	78	79	80

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
1	A1	Uterine cervix	Squamous cell carcinoma	2	F	65
2	A2	Uterine cervix	Squamous cell carcinoma	3	F	57
3	A3	Uterine cervix	Squamous cell carcinoma	2	F	32
4	A4	Uterine cervix	Squamous cell carcinoma	2	F	60
5	A5	Uterine cervix	Squamous cell carcinoma	3	F	41
6	A6	Uterine cervix	Squamous cell carcinoma	3	F	59
7	A7	Uterine cervix	Squamous cell carcinoma	2	F	46
8	A8	Uterine cervix	Squamous cell carcinoma	3	F	34
9	A9	Uterine cervix	Squamous cell carcinoma	2	F	46
10	A10	Uterine cervix	Squamous cell carcinoma	3	F	40
11	B1	Uterine cervix	Squamous cell carcinoma	2	F	49
12	B2	Uterine cervix	Squamous cell carcinoma	2	F	46
13	B3	Uterine cervix	Squamous cell carcinoma	2	F	65
14	B4	Uterine cervix	Squamous cell carcinoma	2	F	53
15	B5	Uterine cervix	Squamous cell carcinoma	2	F	64
16	B6	Uterine cervix	Squamous cell carcinoma	3	F	45
17	B7	Uterine cervix	Squamous cell carcinoma	3	F	53
18	B8	Uterine cervix	Squamous cell carcinoma	2	F	57
19	B9	Uterine cervix	Squamous cell carcinoma	2	F	45
20	B10	Uterine cervix	Squamous cell carcinoma	3	F	39
21	C1	Uterine cervix	Squamous cell carcinoma	3	F	55
22	C2	Uterine cervix	Squamous cell carcinoma	2	F	62



PROTEIN BIOTECHNOLOGIES

23	C3	Uterine cervix	Squamous cell carcinoma	3	F	39
24	C4	Uterine cervix	Squamous cell carcinoma	2	F	53
25	C5	Uterine cervix	Squamous cell carcinoma	3	F	46
26	C6	Uterine cervix	Squamous cell carcinoma	3	F	50
27	C7	Uterine cervix	Squamous cell carcinoma	2	F	35
28	C8	Uterine cervix	Squamous cell carcinoma	3	F	41
29	C9	Uterine cervix	Squamous cell carcinoma	2	F	71
30	C10	Uterine cervix	Squamous cell carcinoma	2	F	35
31	D1	Uterine cervix	Squamous cell carcinoma	1	F	56
32	D2	Uterine cervix	Squamous cell carcinoma	2	F	40
33	D3	Uterine cervix	Squamous cell carcinoma	3	F	71
34	D4	Uterine cervix	Squamous cell carcinoma	3	F	67
35	D5	Uterine cervix	Squamous cell carcinoma	3	F	57
36	D6	Uterine cervix	Squamous cell carcinoma	3	F	31
37	D7	Uterine cervix	Squamous cell carcinoma	2	F	47
38	D8	Uterine cervix	Squamous cell carcinoma	1	F	44
39	D9	Uterine cervix	Squamous cell carcinoma	2	F	63
40	D10	Uterine cervix	Squamous cell carcinoma	2	F	41
41	E1	Uterine cervix	Normal uterine cervix tissue of No. 01	-	F	65
42	E2	Uterine cervix	Normal uterine cervix tissue of No. 02	-	F	57
43	E3	Uterine cervix	Normal uterine cervix tissue of No. 03	-	F	32
44	E4	Uterine cervix	Normal uterine cervix tissue of No. 04	-	F	60
45	E5	Uterine cervix	Normal uterine cervix tissue of No. 05	-	F	41
46	E6	Uterine cervix	Normal uterine cervix tissue of No. 06	-	F	59
47	E7	Uterine cervix	Normal uterine cervix tissue of No. 07	-	F	46
48	E8	Uterine cervix	Normal uterine cervix tissue of No. 08	-	F	34
49	E9	Uterine cervix	Uterine cervix, local squamous cell carcinoma of No. 09	2	F	46
50	E10	Uterine cervix	Cervical canal tissue of No. 10	-	F	40
51	F1	Uterine cervix	Normal uterine cervix tissue	-	F	47
52	F2	Uterine cervix	Normal uterine cervix tissue	-	F	32
53	F3	Uterine cervix	Normal uterine cervix tissue	-	F	44
54	F4	Uterine cervix	Normal uterine cervix tissue	-	F	60
55	F5	Uterine cervix	Normal uterine cervix tissue	-	F	55
56	F6	Uterine cervix	Normal uterine cervix tissue	-	F	26
57	F7	Uterine cervix	Normal uterine cervix tissue	-	F	63
58	F8	Uterine cervix	Normal uterine cervix tissue	-	F	70
59	F9	Uterine cervix	Normal uterine cervix tissue	-	F	41
60	F10	Uterine cervix	Cervical canal tissue	-	F	69



PROTEIN
BIOTECHNOLOGIES

61	G1	Uterine cervix	Normal uterine cervix tissue	-	F	45
62	G2	Uterine cervix	Normal uterine cervix tissue	-	F	43
63	G3	Uterine cervix	Normal uterine cervix tissue	-	F	40
64	G4	Uterine cervix	Normal uterine cervix tissue	-	F	28
65	G5	Uterine cervix	Normal uterine cervix tissue	-	F	39
66	G6	Uterine cervix	Normal uterine cervix tissue	-	F	50
67	G7	Uterine cervix	Cervical canal tissue	-	F	34
68	G8	Uterine cervix	Cervical canal tissue	-	F	42
69	G9	Uterine cervix	Cervical canal tissue with focal squamous metaplasia	-	F	45
70	G10	Uterine cervix	Cervical canal tissue	-	F	48
71	H1	Uterine cervix	Normal uterine cervix tissue	-	F	43
72	H2	Uterine cervix	Cervical canal tissue	-	F	37
73	H3	Uterine cervix	Normal uterine cervix tissue	-	F	59
74	H4	Uterine cervix	Cervical mucosa with focal hyperplasia of parabasal cells	-	F	49
75	H5	Uterine cervix	Normal uterine cervix tissue	-	F	51
76	H6	Uterine cervix	Normal uterine cervix tissue	-	F	31
77	H7	Uterine cervix	Normal uterine cervix tissue	-	F	44
78	H8	Uterine cervix	Normal uterine cervix tissue	-	F	46
79	H9	Uterine cervix	Cervical canal tissue	-	F	49
80	H10	Uterus	Normal uterine cervix tissue	-	F	40

For Research Use Only