



PROTEIN
BIOTECHNOLOGIES

HUMAN THYROID CANCER TISSUE MICROARRAY

Catalog Number: TMA-011

Each core from each tissue represents one specimen selected and pathologically confirmed.
Normal matched or unmatched thyroid tissue is provided for comparison to the cancer tissue.

Cases: 68 **Cores:** 80 **Diameter:** 1.5 mm **Thickness:** 5 µm **Standard IHC:** Thyroglobulin confirmed

Layout: 8 x 10

	1	2	3	4	5	6	7	8	9	10
A	1	2	3	4	5	6	7	8	9	10
B	11	12	13	14	15	16	17	18	19	20
C	21	22	23	24	25	26	27	28	29	30
D	31	32	33	34	35	36	37	38	39	40
E	41	42	43	44	45	46	47	48	49	50
F	51	52	53	54	55	56	57	58	59	60
G	61	62	63	64	65	66	67	68	69	70
H	71	72	73	74	75	76	77	78	79	80

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
1	A1	Thyroid	Follicular carcinoma	-	M	56
2	A2	Thyroid	Papillary carcinoma	-	F	40
3	A3	Thyroid	Papillary carcinoma	-	F	40
4	A4	Thyroid	Papillary carcinoma	-	F	40
5	A5	Thyroid	Papillary carcinoma	-	F	50
6	A6	Thyroid	Follicular carcinoma	-	F	44
7	A7	Thyroid	Follicular carcinoma	-	F	48
8	A8	Thyroid	Papillary carcinoma	-	M	25
9	A9	Thyroid	Papillary carcinoma	-	F	55
10	A10	Thyroid	Follicular papillary carcinoma	-	F	32
11	B1	Thyroid	Papillary carcinoma	-	F	25
12	B2	Thyroid	Papillary carcinoma	-	F	27
13	B3	Thyroid	Papillary carcinoma	-	M	26
14	B4	Thyroid	Papillary carcinoma	-	F	56
15	B5	Thyroid	Papillary carcinoma	-	F	30
16	B6	Thyroid	Papillary carcinoma	-	F	23
17	B7	Thyroid	Papillary carcinoma	-	F	65
18	B8	Thyroid	Papillary carcinoma	-	M	25
19	B9	Thyroid	Papillary carcinoma	-	F	23
20	B10	Thyroid	Papillary carcinoma	-	F	31
21	C1	Thyroid	Papillary carcinoma	-	F	40
22	C2	Thyroid	Papillary carcinoma	-	F	52



PROTEIN

BIOTECHNOLOGIES

23	C3	Thyroid	Papillary carcinoma	-	F	17
24	C4	Thyroid	Papillary carcinoma	-	F	38
25	C5	Thyroid	Papillary carcinoma	-	M	67
26	C6	Thyroid	Papillary carcinoma	-	F	28
27	C7	Thyroid	Papillary carcinoma	-	F	29
28	C8	Thyroid	Papillary carcinoma	-	F	76
29	C9	Thyroid	Papillary carcinoma	-	F	67
30	C10	Thyroid	Papillary carcinoma	-	F	39
31	D1	Thyroid	Papillary carcinoma	-	M	36
32	D2	Thyroid	Follicular papillary carcinoma	-	F	21
33	D3	Thyroid	Follicular papillary carcinoma	-	F	23
34	D4	Thyroid	Follicular papillary carcinoma	-	F	65
35	D5	Thyroid	Follicular papillary carcinoma	-	F	26
36	D6	Thyroid	Sarcomatoid carcinoma	-	F	40
37	D7	Thyroid	Squamous cell carcinoma	3	F	56
38	D8	Thyroid	Medullary carcinoma	-	M	49
39	D9	Thyroid	Medullary carcinoma	-	F	39
40	D10	Thyroid	Malignant lymphoma	-	F	63
41	E1	Thyroid	Normal thyroid tissue of No.01	-	M	56
42	E2	Thyroid	Normal thyroid tissue of No.02	-	F	40
43	E3	Thyroid	Focal lymphocyte infiltrating in tissue No.03	-	F	40
44	E4	Thyroid	Normal thyroid tissue of No.04	-	F	40
45	E5	Thyroid	Normal thyroid tissue of No.05	-	F	50
46	E6	Thyroid	Normal thyroid tissue of No.06	-	F	44
47	E7	Thyroid	Normal thyroid tissue of No.07	-	F	48
48	E8	Thyroid	Normal thyroid tissue of No.08	-	M	25
49	E9	Thyroid	Normal thyroid tissue of No.09	-	F	55
50	E10	Thyroid	Normal thyroid tissue of No.10	-	F	32
51	F1	Thyroid	Normal thyroid tissue of No.11	-	F	25
52	F2	Thyroid	Normal thyroid tissue of No.12	-	F	27
53	F3	Thyroid	Normal thyroid tissue	-	F	48
54	F4	Thyroid	Normal thyroid tissue	-	F	42
55	F5	Thyroid	Normal thyroid tissue	-	F	36
56	F6	Thyroid	Normal thyroid tissue	-	F	39
57	F7	Thyroid	Normal thyroid tissue	-	F	36
58	F8	Thyroid	Normal thyroid tissue	-	F	30
59	F9	Thyroid	Normal thyroid tissue	-	M	48
60	F10	Thyroid	Normal thyroid tissue	-	M	39



PROTEIN

BIOTECHNOLOGIES

61	G1	Thyroid	Normal thyroid tissue	-	F	31
62	G2	Thyroid	Normal thyroid tissue	-	F	51
63	G3	Thyroid	Normal thyroid tissue	-	F	46
64	G4	Thyroid	Normal thyroid tissue	-	F	44
65	G5	Thyroid	Normal thyroid tissue	-	F	57
66	G6	Thyroid	Normal thyroid tissue	-	F	52
67	G7	Thyroid	Normal thyroid tissue	-	F	20
68	G8	Thyroid	Normal thyroid tissue	-	F	37
69	G9	Thyroid	Normal thyroid tissue	-	F	37
70	G10	Thyroid	Normal thyroid tissue	-	F	32
71	H1	Thyroid	Normal thyroid tissue	-	M	27
72	H2	Thyroid	Normal thyroid tissue	-	F	43
73	H3	Thyroid	Normal thyroid tissue	-	F	38
74	H4	Thyroid	Normal thyroid tissue	-	F	62
75	H5	Thyroid	Normal thyroid tissue	-	F	43
76	H6	Thyroid	Normal thyroid tissue	-	M	55
77	H7	Thyroid	Normal thyroid tissue	-	F	40
78	H8	Thyroid	Normal thyroid tissue	-	F	65
79	H9	Thyroid	Normal thyroid tissue	-	F	31
80	H10	Thyroid	Normal thyroid tissue	-	F	57

For Research Use Only