



PROTEIN
BIOTECHNOLOGIES

HUMAN LARYNGEAL & PHARYNGEAL CANCER TISSUE MICROARRAY

Catalog Number: TMA-017

Each core from each tissue represents one specimen selected and pathologically confirmed. Normal matched or unmatched tissue is provided for comparison to the cancer tissue.

Cases: 40 **Cores:** 80 **Diameter:** 1.5 mm **Thickness:** 5 µm **Standard IHC:** Cytokeratin confirmed

Layout: 8 x 10

	1	2	3	4	5	6	7	8	9	10
A	1	2	3	4	5	6	7	8	9	10
B	11	12	13	14	15	16	17	18	19	20
C	21	22	23	24	25	26	27	28	29	30
D	31	32	33	34	35	36	37	38	39	40
E	41	42	43	44	45	46	47	48	49	50
F	51	52	53	54	55	56	57	58	59	60
G	61	62	63	64	65	66	67	68	69	70
H	71	72	73	74	75	76	77	78	79	80

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
1	A1	Larynx	Squamous cell carcinoma	3	M	56
2	A2	Larynx	Squamous cell carcinoma	1	M	51
3	A3	Larynx	Squamous cell carcinoma	3	M	53
4	A4	Larynx	Squamous cell carcinoma	2	M	68
5	A5	Larynx	Squamous cell carcinoma	2	M	56
6	A6	Larynx	Squamous cell carcinoma	3	M	58
7	A7	Larynx	Squamous cell carcinoma	1	M	47
8	A8	Larynx	Adenosquamous carcinoma	-	M	73
9	A9	Larynx	Squamous cell carcinoma	1	M	74
10	A10	Larynx	Squamous cell carcinoma	2	M	65
11	B1	Nose	Papillary adenocarcinoma	1	F	15
12	B2	Larynx	Squamous cell carcinoma	1	M	58
13	B3	Larynx	Squamous cell carcinoma	1	M	61
14	B4	Larynx	Squamous cell carcinoma	1	M	65
15	B5	Larynx	Squamous cell carcinoma	1	M	48
16	B6	Larynx	Squamous cell carcinoma	2	M	50
17	B7	Larynx	Squamous cell carcinoma	2	M	45
18	B8	Larynx	Squamous cell carcinoma	3	M	63
19	B9	Larynx	Squamous cell carcinoma	1	M	65
20	B10	Larynx	Squamous cell carcinoma	1	M	72
21	C1	Larynx	Squamous cell carcinoma	2	M	71
22	C2	Larynx	Squamous cell carcinoma	1	M	47



PROTEIN BIOTECHNOLOGIES

23	C3	Larynx	Squamous cell carcinoma	1	M	45
24	C4	Larynx	Squamous cell carcinoma	1	M	56
25	C5	Larynx	Squamous cell carcinoma	2	M	58
26	C6	Larynx	Basaloid squamous cell carcinoma	-	M	58
27	C7	Larynx	Squamous cell carcinoma	1	M	56
28	C8	Larynx	Squamous cell carcinoma	2	M	67
29	C9	Larynx	Squamous cell carcinoma	1	M	54
30	C10	Larynx	Squamous cell carcinoma	1	M	60
31	D1	Larynx	Squamous cell carcinoma	1	M	65
32	D2	Nose	Squamous cell carcinoma	1	M	33
33	D3	Larynx	Squamous cell carcinoma	2	M	62
34	D4	Larynx	Squamous cell carcinoma	1	M	65
35	D5	Larynx	Squamous cell carcinoma	1	F	69
36	D6	Larynx	Squamous cell carcinoma	2	M	59
37	D7	Larynx	Salivary gland tissue	-	M	68
38	D8	Larynx	Squamous cell carcinoma	3	M	67
39	D9	Larynx	Squamous cell carcinoma	2	M	41
40	D10	Larynx	Squamous cell carcinoma	2	M	72
41	E1	Larynx	Squamous cell carcinoma	1	M	62
42	E2	Larynx	Squamous cell carcinoma	1	M	63
43	E3	Larynx	Squamous cell carcinoma	3	M	70
44	E4	Larynx	Squamous cell carcinoma	3	M	64
45	E5	Larynx	Squamous cell carcinoma	3	M	67
46	E6	Larynx	Squamous cell carcinoma	2	M	54
47	E7	Larynx	Squamous cell carcinoma	3	M	48
48	E8	Larynx	Squamous cell carcinoma	3	M	68
49	E9	Larynx	Squamous cell carcinoma	2	M	63
50	E10	Larynx	Squamous cell carcinoma	3	M	68
51	F1	Larynx	Squamous cell carcinoma	2	M	56
52	F2	Larynx	Squamous cell carcinoma	3	M	60
53	F3	Larynx	Squamous cell carcinoma	3	M	47
54	F4	Larynx	Squamous cell carcinoma	2	F	50
55	F5	Larynx	Squamous cell carcinoma	2	M	58
56	F6	Larynx	Squamous cell carcinoma	2	F	73
57	F7	Larynx	Squamous cell carcinoma	3	M	51
58	F8	Larynx	Squamous cell carcinoma	2	M	47
59	F9	Larynx	Squamous cell carcinoma	2	M	70
60	F10	Larynx	Squamous cell carcinoma	3	M	52



PROTEIN
BIOTECHNOLOGIES

61	G1	Pharynx	Diffuse B cell lymphoma	-	M	60
62	G2	Nose	Diffuse B cell lymphoma	-	M	74
63	G3	Tongue	T cell lymphoma	-	F	73
64	G4	Larynx	Basaloid squamous cell carcinoma	-	M	63
65	G5	Larynx	Squamous cell carcinoma	2	M	74
66	G6	Tongue	Squamous cell carcinoma	1	M	57
67	G7	Nose	Squamous cell carcinoma	3	F	56
68	G8	Larynx	Fibrosarcoma	-	F	39
69	G9	Tongue	Acinic cell carcinoma	-	F	75
70	G10	Tongue	Squamous cell carcinoma	2	F	24
71	H1	Larynx	Normal larynx tissue of No. 01	-	M	56
72	H2	Larynx	Normal larynx tissue of No. 02	-	M	51
73	H3	Larynx	Normal larynx tissue of No. 03	-	M	53
74	H4	Larynx	Normal larynx tissue of No. 04	-	M	68
75	H5	Larynx	Normal larynx tissue of No. 05	-	M	56
76	H6	Larynx	Normal larynx tissue of No. 06	-	M	58
77	H7	Larynx	Chronic mucosal inflammation of No. 07	-	M	47
78	H8	Larynx	Normal larynx tissue of No. 08	-	M	73
79	H9	Larynx	Normal larynx tissue of No. 09	-	M	74
80	H10	Larynx	Normal larynx tissue of No. 10	-	M	65

For Research Use Only