



PROTEIN
BIOTECHNOLOGIES

HUMAN LUNG CANCER TISSUE ARRAY

Catalog Number: TMA-031

Each three (3) dots from three different tissue spots represents one single specimen that was selected and pathologically confirmed (Lung squamous cell carcinoma tissue array)

Cases: 21 **Cores:** 63 **Diameter:** 1.5 mm **Thickness:** 5 µm **Standard IHC:** H-CK confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Lung	Squamous cell carcinoma	II	M	58
02	A2	Lung	Squamous cell carcinoma	II	M	58
03	A3	Lung	Squamous cell carcinoma	II	M	58
04	A4	Lung	Squamous cell carcinoma	II	M	55
05	A5	Lung	Squamous cell carcinoma	II	M	55
06	A6	Lung	Squamous cell carcinoma	II	M	55
07	A7	Lung	Squamous cell carcinoma	II	M	46
08	A8	Lung	Squamous cell carcinoma	II	M	46
09	A9	Lung	Squamous cell carcinoma	II	M	46
10	B1	Lung	Normal lung tissue		M	65
11	B2	Lung	Normal lung tissue		M	65
12	B3	Lung	Normal lung tissue		M	65
13	B4	Lung	Squamous cell carcinoma	II	M	55
14	B5	Lung	Squamous cell carcinoma	II	M	55
15	B6	Lung	Squamous cell carcinoma	II	M	55
16	B7	Lung	Few carcinoma cells		M	68
17	B8	Lung	Squamous cell carcinoma	II	M	68
18	B9	Lung	Squamous cell carcinoma	II	M	68
19	C1	Lung	Squamous cell carcinoma	II	M	69
20	C2	Lung	Squamous cell carcinoma	II	M	69
21	C3	Lung	Squamous cell carcinoma	II	M	69
22	C4	Lung	Squamous cell carcinoma	II	M	49
23	C5	Lung	Squamous cell carcinoma	II	M	49
24	C6	Lung	Squamous cell carcinoma	II	M	49
25	C7	Lung	Squamous cell carcinoma	II	M	62
26	C8	Lung	Squamous cell carcinoma	II	M	62
27	C9	Lung	Squamous cell carcinoma	II	M	62
28	D1	Lung	Squamous cell carcinoma	II	M	64
29	D2	Lung	Squamous cell carcinoma	II	M	64
30	D3	Lung	Squamous cell carcinoma	II	M	64
31	D4	Lung	Squamous cell carcinoma	II	M	62
32	D5	Lung	Squamous cell carcinoma	II	M	62
33	D6	Lung	Squamous cell carcinoma	II	M	62
34	D7	Lung	Squamous cell carcinoma	II	M	73



PROTEIN BIOTECHNOLOGIES

35	D8	Lung	Squamous cell carcinoma	II	M	73
36	D9	Lung	Squamous cell carcinoma	II	M	73
37	E1	Lung	Squamous cell carcinoma	II	M	42
38	E2	Lung	Squamous cell carcinoma	II	M	42
39	E3	Lung	Squamous cell carcinoma	II	M	42
40	E4	Lung	Squamous cell carcinoma	II	F	51
41	E5	Lung	Squamous cell carcinoma	II	F	51
42	E6	Lung	Squamous cell carcinoma	II	F	51
43	E7	Lung	Squamous cell carcinoma	II	M	35
44	E8	Lung	Squamous cell carcinoma	II	M	35
45	E9	Lung	Squamous cell carcinoma	II	M	35
46	F1	Lung	Squamous cell carcinoma	II	M	45
47	F2	Lung	Squamous cell carcinoma	II	M	45
48	F3	Lung	Squamous cell carcinoma	II	M	45
49	F4	Lung	Squamous cell carcinoma	II	M	67
50	F5	Lung	Squamous cell carcinoma	II	M	67
51	F6	Lung	Squamous cell carcinoma	II	M	67
52	F7	Lung	Yellow tumor-like lesion of margin of carcinoma		M	70
53	F8	Lung	Yellow tumor-like lesion of margin of carcinoma		M	70
54	F9	Lung	Yellow tumor-like lesion of margin of carcinoma		M	70
55	G1	Lung	Squamous cell carcinoma	II	M	59
56	G2	Lung	Squamous cell carcinoma	II	M	59
57	G3	Lung	Squamous cell carcinoma	II	M	59
58	G4	Lung	Squamous cell carcinoma	II	M	69
59	G5	Lung	Squamous cell carcinoma	II	M	69
60	G6	Lung	Squamous cell carcinoma	II	M	69
61	G7	Lung	Normal lung tissue		M	70
62	G8	Lung	Normal lung tissue		M	70
63	G9	Lung	Normal lung tissue		M	70

For Research Use Only