



PROTEIN
BIOTECHNOLOGIES

HUMAN LUNG CANCER TISSUE ARRAY

Catalog Number: TMA-035

Each two (2) dots from two different tissue spots represents one single specimen that was selected and pathologically confirmed (Lung small cell carcinoma tissue array)

Cases: 33 **Cores:** 63 **Diameter:** 1.5 mm **Thickness:** 5 µm **Standard IHC:** L-CK confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Lung	Small cell carcinoma		F	69
02	A2	Lung	Small cell carcinoma		F	69
03	A3	Lung	Small cell carcinoma		F	42
04	A4	Lung	Small cell carcinoma		F	42
05	A5	Lung	Small cell carcinoma		M	69
06	A6	Lung	Small cell carcinoma		M	69
07	A7	Lung	Small cell carcinoma		M	57
08	A8	Lung	Small cell carcinoma		M	57
09	A9	Lung	Small cell carcinoma		M	70
10	B1	Lung	Small cell carcinoma		M	61
11	B2	Lung	Small cell carcinoma		M	61
12	B3	Lung	Small cell carcinoma		M	50
13	B4	Lung	Small cell carcinoma		M	50
14	B5	Lung	Small cell carcinoma		M	62
15	B6	Lung	Small cell carcinoma		M	62
16	B7	Lung	Small cell carcinoma		F	50
17	B8	Lung	Small cell carcinoma		F	50
18	B9	Lung	Small cell carcinoma		M	70
19	C1	Lung	Small cell carcinoma		M	60
20	C2	Lung	Small cell carcinoma		M	60
21	C3	Lung	Small cell carcinoma		F	53
22	C4	Lung	Small cell carcinoma		F	53
23	C5	Lung	Small cell carcinoma		M	57
24	C6	Lung	Small cell carcinoma		M	57
25	C7	Lung	Small cell carcinoma		M	64
26	C8	Lung	Small cell carcinoma		M	64
27	C9	Lung	Small cell carcinoma		F	51
28	D1	Lung	Small cell carcinoma		F	46
29	D2	Lung	Small cell carcinoma		F	46
30	D3	Lung	Small cell carcinoma		M	54
31	D4	Lung	Small cell carcinoma		M	54
32	D5	Lung	Small cell carcinoma		M	72
33	D6	Lung	Small cell carcinoma		M	72
34	D7	Lung	Small cell carcinoma		M	73
35	D8	Lung	Small cell carcinoma		M	73



PROTEIN

BIOTECHNOLOGIES

36	D9	Lung	Small cell carcinoma		F	51
37	E1	Lung	Small cell carcinoma		F	66
38	E2	Lung	Small cell carcinoma		F	66
39	E3	Lung	Small cell carcinoma		F	56
40	E4	Lung	Small cell carcinoma		F	56
41	E5	Lung	Small cell carcinoma		F	33
42	E6	Lung	Small cell carcinoma		F	33
43	E7	Lung	Small cell carcinoma		M	29
44	E8	Lung	Small cell carcinoma		M	29
45	E9	Lung	Small cell carcinoma		F	65
46	F1	Lung	Small cell carcinoma		M	53
47	F2	Lung	Small cell carcinoma		M	53
48	F3	Lung	Small cell carcinoma		F	62
49	F4	Lung	Small cell carcinoma		F	62
50	F5	Lung	Small cell carcinoma		M	55
51	F6	Lung	Small cell carcinoma		M	55
52	F7	Lung	Small cell carcinoma		F	56
53	F8	Lung	Small cell carcinoma		F	56
54	F9	Lung	Small cell carcinoma		F	65
55	G1	Lung	Small cell carcinoma		M	40
56	G2	Lung	Small cell carcinoma		M	40
57	G3	Lung	Small cell carcinoma		M	55
58	G4	Lung	Small cell carcinoma		M	55
59	G5	Lung	Small cell carcinoma		F	65
60	G6	Lung	Small cell carcinoma		F	65
61	G7	Lung	Normal lung tissue		M	65
62	G8	Lung	Normal lung tissue		M	42
63	G9	Lung	Normal lung tissue		M	11

For Research Use Only