



PROTEIN
BIOTECHNOLOGIES

HUMAN COLON CANCER TISSUE ARRAY

Catalog Number: TMA-054

Each three (3) dots from three different tissue spots represents one single specimen selected and pathologically confirmed (normal colon tissue array).

Cases: 23 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** Actin confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Colon	Normal colon mucosae tissue	-	F	56
02	A2	Colon	Normal colon mucosae tissue	-	F	56
03	A3	Colon	Fibred tissue	-	F	56
04	A4	Colon	Normal colon mucosae tissue	-	M	51
05	A5	Colon	Normal colon mucosae tissue	-	M	51
06	A6	Colon	Normal colon mucosae tissue	-	M	51
07	A7	Colon	Normal colon mucosae tissue	-	F	50
08	A8	Colon	Normal colon mucosae tissue	-	F	50
09	A9	Colon	Smooth muscle	-	F	50
10	B1	Colon	Normal colon mucosae tissue	-	F	51
11	B2	Colon	Normal colon mucosae tissue	-	F	51
12	B3	Colon	Normal colon mucosae tissue	-	F	51
13	B4	Colon	Normal colon mucosae tissue	-	M	70
14	B5	Colon	Normal colon mucosae tissue	-	M	70
15	B6	Colon	Normal colon mucosae tissue	-	M	70
16	B7	Colon	Normal colon mucosae tissue	-	M	26
17	B8	Colon	Normal colon mucosae tissue	-	M	26
18	B9	Colon	Normal colon mucosae tissue	-	M	26
19	C1	Colon	Smooth muscle	-	M	50
20	C2	Colon	Smooth muscle	-	M	50
21	C3	Colon	Smooth muscle	-	M	50
22	C4	Colon	Smooth muscle	-	M	49
23	C5	Colon	Normal colon mucosae tissue	-	M	49
24	C6	Colon	Normal colon mucosae tissue	-	M	49
25	C7	Colon	Normal colon mucosae tissue	-	M	27
26	C8	Colon	Normal colon mucosae tissue	-	M	27
27	C9	Colon	Fibred tissue	-	M	27
28	D1	Colon	Normal colon mucosae tissue	-	M	38
29	D2	Colon	Normal colon mucosae tissue	-	M	38
30	D3	Colon	Normal colon mucosae tissue	-	M	38
31	D4	Colon	Normal colon mucosae tissue	-	M	68
32	D5	Colon	Normal colon mucosae tissue	-	M	68
33	D6	Colon	Normal colon mucosae tissue	-	M	68
34	D7	Colon	Fibred fat tissue	-	M	38



PROTEIN

BIOTECHNOLOGIES

35	D8	Colon	Normal colon mucosae tissue	-	M	38
36	D9	Colon	Normal colon mucosae tissue	-	M	38
37	E1	Colon	Normal colon mucosae tissue	-	F	55
38	E2	Colon	Normal colon mucosae tissue	-	F	55
39	E3	Colon	Normal colon mucosae tissue	-	F	55
40	E4	Colon	Normal colon mucosae tissue	-	M	61
41	E5	Colon	Smooth muscle fibred	-	M	61
42	E6	Colon	Smooth muscle	-	M	61
43	E7	Colon	Normal colon mucosae tissue	-	F	52
44	E8	Colon	Normal colon mucosae tissue	-	F	52
45	E9	Colon	Normal colon mucosae tissue	-	F	52
46	F1	Colon	Fibred fat tissue	-	F	11
47	F2	Colon	Fibred fat tissue	-	F	11
48	F3	Colon	Fibred fat tissue	-	F	11
49	F4	Colon	Normal colon mucosae tissue	-	M	64
50	F5	Colon	Fibred tissue	-	M	64
51	F6	Colon	Normal colon mucosae tissue	-	M	64
52	F7	Colon	Normal colon mucosae tissue	-	M	65
53	F8	Colon	Normal colon mucosae tissue	-	M	65
54	F9	Colon	Normal colon mucosae tissue	-	M	65
55	G1	Colon	Normal colon mucosae tissue	-	M	69
56	G2	Colon	Smooth muscle	-	M	69
57	G3	Colon	Normal tissue	-	M	69
58	G4	Colon	Normal colon mucosae tissue	-	M	82
59	G5	Colon	Smooth muscle	-	M	82
60	G6	Colon	Normal colon mucosae tissue	-	M	82
61	G7	Colon	Mucinous adenocarcinoma	-	F	51
62	G8	Colon	Adenocarcinoma	-	M	50
63	G9	Colon	Adenocarcinoma	-	M	82

For Research Use Only