



PROTEIN
BIOTECHNOLOGIES

HUMAN COLON CANCER TISSUE ARRAY

Catalog Number: TMA-056

Each three (3) dots from three different tissue spots represents one single specimen selected and pathologically confirmed.(Colon Adenocarcinoma tissue array)

Cases: 23 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** P53 confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Colon	Adenocarcinoma	III	M	37
02	A2	Colon	Adenocarcinoma	III	M	37
03	A3	Colon	Adenocarcinoma	I	M	37
04	A4	Colon	Adenocarcinoma	I	M	70
05	A5	Colon	Adenocarcinoma	I	M	70
06	A6	Colon	Adenocarcinoma	I	M	70
07	A7	Colon	Adenocarcinoma	I	F	62
08	A8	Colon	Adenocarcinoma	I	F	62
09	A9	Colon	Adenocarcinoma	I	F	62
10	B1	Colon	Adenocarcinoma	I	M	49
11	B2	Colon	Adenocarcinoma	I	M	49
12	B3	Colon	Adenocarcinoma	I	M	49
13	B4	Colon	Adenocarcinoma	I	M	70
14	B5	Colon	Adenocarcinoma	I	M	70
15	B6	Colon	Adenocarcinoma	I	M	70
16	B7	Colon	Adenocarcinoma	I	M	82
17	B8	Colon	Adenocarcinoma	I	M	82
18	B9	Colon	Adenocarcinoma	II	M	82
19	C1	Colon	Adenocarcinoma	I	M	38
20	C2	Colon	Adenocarcinoma	I	M	38
21	C3	Colon	Adenocarcinoma	I	M	38
22	C4	Colon	Adenocarcinoma	II	F	50
23	C5	Colon	Adenocarcinoma	II	F	50
24	C6	Colon	Adenocarcinoma	III	F	50
25	C7	Colon	Adenocarcinoma	III	M	50
26	C8	Colon	Adenocarcinoma	II	M	50
27	C9	Colon	Adenocarcinoma	II	M	50
28	D1	Colon	Adenocarcinoma	I	M	57
29	D2	Colon	Adenocarcinoma	I	M	57
30	D3	Colon	Adenocarcinoma	I	M	57
31	D4	Colon	Adenocarcinoma	II	F	72
32	D5	Colon	Adenocarcinoma	II	F	72
33	D6	Colon	Adenocarcinoma	III	F	72
34	D7	Colon	Adenocarcinoma	I	M	26
35	D8	Colon	Adenocarcinoma	I	M	26



PROTEIN

BIOTECHNOLOGIES

36	D9	Colon	Adenocarcinoma	I	M	26
37	E1	Colon	Inflammatory tissue	-	F	55
38	E2	Colon	Mucinous adenocarcinoma	I	F	55
39	E3	Colon	Mucinous adenocarcinoma	I	F	55
40	E4	Colon	Adenocarcinoma	III	F	48
41	E5	Colon	Adenocarcinoma	III	F	48
42	E6	Colon	Adenocarcinoma	III	F	48
43	E7	Colon	Adenocarcinoma	I	M	86
44	E8	Colon	Adenocarcinoma	I	M	86
45	E9	Colon	Adenocarcinoma	I	M	86
46	F1	Colon	Inflammatory tissue	-	M	64
47	F2	Colon	Mucinous adenocarcinoma	-	M	64
48	F3	Colon	Mucinous adenocarcinoma	-	M	64
49	F4	Colon	Adenocarcinoma	I	M	77
50	F5	Colon	Adenocarcinoma	I	M	77
51	F6	Colon	Adenocarcinoma	I	M	77
52	F7	Colon	Adenocarcinoma	I	M	62
53	F8	Colon	Adenocarcinoma	I	M	62
54	F9	Colon	Adenocarcinoma	I	M	62
55	G1	Colon	Adenocarcinoma	I	F	33
56	G2	Colon	Adenocarcinoma	I	F	33
57	G3	Colon	Adenocarcinoma	I	F	33
58	G4	Colon	Adenocarcinoma	I	M	66
59	G5	Colon	Adenocarcinoma	I	M	43
60	G6	Colon	Adenocarcinoma	I	M	52
61	G7	Colon	Normal colon mucosae tissue	-	F	52
62	G8	Colon	Normal colon mucosae tissue	-	F	11
63	G9	Colon	Normal colon mucosae tissue	-	F	56

For Research Use Only