



PROTEIN
BIOTECHNOLOGIES

HUMAN COLON CANCER TISSUE ARRAY

Catalog Number: TMA-057

Each three (3) dots from three different tissue spots represents one single specimen selected and pathologically confirmed.(Colon Adenocarcinoma tissue array)

Cases: 23 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** L-CK confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Colon	Adenocarcinoma	I	F	62
02	A2	Colon	Adenocarcinoma	I	F	62
03	A3	Colon	Adenocarcinoma	I	F	62
04	A4	Colon	Adenocarcinoma	I	M	26
05	A5	Colon	Adenocarcinoma	I	M	26
06	A6	Colon	Adenocarcinoma	I	M	26
07	A7	Colon	Mucous adenocarcinoma	I	F	51
08	A8	Colon	Adenocarcinoma	I	F	51
09	A9	Colon	Mucous adenocarcinoma	I	F	51
10	B1	Colon	Adenocarcinoma	I	M	55
11	B2	Colon	Adenocarcinoma	I	M	55
12	B3	Colon	Adenocarcinoma	I	M	55
13	B4	Colon	Adenocarcinoma	I	M	59
14	B5	Colon	Adenocarcinoma	I	M	59
15	B6	Colon	Adenocarcinoma	I	M	59
16	B7	Colon	Adenocarcinoma	I	F	57
17	B8	Colon	Adenocarcinoma	I	F	57
18	B9	Colon	Adenocarcinoma	I	F	57
19	C1	Colon	Adenocarcinoma	II	M	86
20	C2	Colon	Adenocarcinoma	II	M	86
21	C3	Colon	Adenocarcinoma	II	M	86
22	C4	Colon	Adenocarcinoma	I	F	53
23	C5	Colon	Adenocarcinoma	I	F	53
24	C6	Colon	Adenocarcinoma	I	F	53
25	C7	Colon	Smooth muscle and a few carcinoma cells infiltrating	-	M	20
26	C8	Colon	Adenocarcinoma	III	M	20
27	C9	Colon	Adenocarcinoma	III	M	20
28	D1	Colon	Adenocarcinoma	I	M	75
29	D2	Colon	Adenocarcinoma	I	M	75
30	D3	Colon	Adenocarcinoma	I	M	75
31	D4	Colon	Adenocarcinoma	I	M	81
32	D5	Colon	Adenocarcinoma	I	M	81
33	D6	Colon	Adenocarcinoma	I	M	81
34	D7	Colon	Adenocarcinoma	I	M	72
35	D8	Colon	Adenocarcinoma	I	M	72



PROTEIN BIOTECHNOLOGIES

36	D9	Colon	Adenocarcinoma	I	M	72
37	E1	Colon	Adenocarcinoma	III	M	48
38	E2	Colon	A few carcinoma cells infiltrating	-	M	48
39	E3	Colon	Adenocarcinoma	III	M	48
40	E4	Colon	Adenocarcinoma	I	F	59
41	E5	Colon	Adenocarcinoma	I	F	59
42	E6	Colon	Adenocarcinoma	I	F	59
43	E7	Colon	Adenocarcinoma	I	M	44
44	E8	Colon	Adenocarcinoma	I	M	44
45	E9	Colon	Adenocarcinoma	I	M	44
46	F1	Colon	Adenocarcinoma	I	F	67
47	F2	Colon	Adenocarcinoma	I	F	67
48	F3	Colon	Adenocarcinoma	I	F	67
49	F4	Colon	Adenocarcinoma	I	F	64
50	F5	Colon	Adenocarcinoma	I	F	64
51	F6	Colon	Adenocarcinoma	I	F	64
52	F7	Colon	Adenocarcinoma	I	M	60
53	F8	Colon	Adenocarcinoma	I	M	60
54	F9	Colon	Adenocarcinoma	I	M	60
55	G1	Colon	Adenocarcinoma	I	M	47
56	G2	Colon	Adenocarcinoma	I	M	47
57	G3	Colon	Adenocarcinoma	I	M	47
58	G4	Colon	Adenocarcinoma	I	M	62
59	G5	Colon	Adenocarcinoma	I	M	62
60	G6	Colon	Adenocarcinoma	I	M	62
61	G7	Intestine	Mucosae of small intestine	-	M	69
62	G8	Colon	Normal colon mucosae tissue	-	F	56
63	G9	Colon	Normal colon mucosae tissue	-	M	70

For Research Use Only