



PROTEIN
BIOTECHNOLOGIES

HUMAN COLON CANCER TISSUE ARRAY

Catalog Number: TMA-058

Each three (3) dots from three different tissue spots represents one single specimen that was selected and pathologically confirmed (Colon Adenocarcinoma tissue array)

Cases: 21 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** L-CK confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Locaton	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Colon	Adenocarcinoma	II	M	37
02	A2	Colon	Adenocarcinoma	II	M	37
03	A3	Colon	Adenocarcinoma	II	M	37
04	A4	Colon	Adenocarcinoma	I	F	52
05	A5	Colon	Adenocarcinoma	I	F	52
06	A6	Colon	Adenocarcinoma	I	F	52
07	A7	Colon	Adenocarcinoma	II	M	61
08	A8	Colon	Adenocarcinoma	I	M	61
09	A9	Colon	Adenocarcinoma	II	M	61
10	B1	Colon	Adenocarcinoma	I	F	48
11	B2	Colon	Adenocarcinoma	I	F	48
12	B3	Colon	Adenocarcinoma	I	F	48
13	B4	Colon	Adenocarcinoma	I	M	82
14	B5	Colon	Adenocarcinoma	I	M	82
15	B6	Colon	Adenocarcinoma	I	M	82
16	B7	Colon	Adenocarcinoma	II	M	75
17	B8	Colon	Adenocarcinoma	I	M	75
18	B9	Colon	Adenocarcinoma	I	M	75
19	C1	Colon	Adenocarcinoma	I	F	43
20	C2	Colon	Adenocarcinoma	I	F	43
21	C3	Colon	Adenocarcinoma	I	F	43
22	C4	Colon	Adenocarcinoma	I	M	55
23	C5	Colon	Adenocarcinoma	I	M	55
24	C6	Colon	Adenocarcinoma	II	M	55
25	C7	Colon	Adenocarcinoma	I	F	51
26	C8	Colon	Adenocarcinoma	I	F	51
27	C9	Colon	Adenocarcinoma	I	F	51
28	D1	Colon	Adenocarcinoma	I	M	60
29	D2	Colon	Adenocarcinoma	I	M	60
30	D3	Colon	Adenocarcinoma	I	M	60
31	D4	Colon	Adenocarcinoma	II	F	17
32	D5	Colon	Adenocarcinoma	II	F	17
33	D6	Colon	Adenocarcinoma	II	F	17
34	D7	Colon	Adenocarcinoma	I	M	62
35	D8	Colon	Adenocarcinoma	I	M	62



PROTEIN

BIOTECHNOLOGIES

36	D9	Colon	Adenocarcinoma	II	M	62
37	E1	Colon	Adenocarcinoma	II	M	54
38	E2	Colon	Adenocarcinoma	II	M	54
39	E3	Colon	Adenocarcinoma	II	M	54
40	E4	Colon	Adenocarcinoma	II	M	66
41	E5	Colon	Adenocarcinoma	II	M	66
42	E6	Colon	Adenocarcinoma	II	M	66
43	E7	Colon	Adenocarcinoma	I	F	41
44	E8	Colon	Adenocarcinoma	I	F	41
45	E9	Colon	Adenocarcinoma	I	F	41
46	F1	Colon	Adenocarcinoma	II	F	63
47	F2	Colon	Adenocarcinoma	II	F	63
48	F3	Colon	Adenocarcinoma	II	F	63
49	F4	Colon	Adenocarcinoma	II	M	68
50	F5	Colon	Adenocarcinoma	II	M	68
51	F6	Colon	Smooth muscle and chronic inflammation	-	M	68
52	F7	Colon	Adenocarcinoma	II	M	46
53	F8	Colon	Adenocarcinoma	II	M	46
54	F9	Colon	Adenocarcinoma	II	M	46
55	G1	Colon	Adenocarcinoma	I	F	68
56	G2	Colon	Adenocarcinoma	I	F	68
57	G3	Colon	Adenocarcinoma	I	F	68
58	G4	Colon	Adenocarcinoma	I	M	68
59	G5	Colon	Adenocarcinoma	I	M	68
60	G6	Colon	Adenocarcinoma	II	M	68
61	G7	Colon	Smooth muscle and fibrous tissues	-	F	52
62	G8	Colon	Mucosae of colon	-	M	61
63	G9	Intestine	Mucosae of small intestine	-	M	50

For Research Use Only