



PROTEIN
BIOTECHNOLOGIES

HUMAN UTERUS CANCER TISSUE ARRAY

Catalog Number: TMA-069

Each dot represents a DISEASED tissue spot from one individual specimen that was selected and pathologically confirmed

Cases: 63 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** CK confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Uterus	Carcinoma tissue	-	F	61
02	A2	Uterus	Endometrial adenocarcinoma	II	F	52
03	A3	Uterus	Endometrial adenocarcinoma	I	F	64
04	A4	Uterus	Endometrial adenocarcinoma	I	F	52
05	A5	Uterus	Endometrial adenocarcinoma	I	F	51
06	A6	Uterus	Endometrial adenocarcinoma	III	F	58
07	A7	Uterus	Endometrial adenocarcinoma	II	F	44
08	A8	Uterus	Endometrial adenocarcinoma	II	F	54
09	A9	Uterus	Endometrial adenocarcinoma	I	F	38
10	B1	Uterus	Endometrial adenocarcinoma	I	F	48
11	B2	Uterus	Invasive hydatidiform mole	-	F	55
12	B3	Uterus	Endometrial adenocarcinoma	II	F	60
13	B4	Uterus	Endometrial adenocarcinoma	II	F	63
14	B5	Uterus	Endometrial adenocarcinoma	II	F	53
15	B6	Uterus	Endometrial adenocarcinoma	I	F	53
16	B7	Uterus	Endometrial adenocarcinoma	II	F	66
17	B8	Uterus	Endometrial adenocarcinoma	II	F	56
18	B9	Uterus	Endometrial adenocarcinoma	II	F	49
19	C1	Uterus	Endometrial adenocarcinoma	II	F	54
20	C2	Uterus	Endometrial adenocarcinoma	I	F	58
21	C3	Uterus	Endometrial adenocarcinoma	I	F	62
22	C4	Uterus	Endometrial adenocarcinoma	II	F	57
23	C5	Uterus	Endometrial adenocarcinoma	II	F	51
24	C6	Uterus	Endometrial adenocarcinoma	II	F	45
25	C7	Uterus	Endometrial adenocarcinoma	I	F	55
26	C8	Uterus	Endometrial adenocarcinoma	II	F	59
27	C9	Uterus	Endometrial adenocarcinoma	II	F	51
28	D1	Uterus	Endometrial adenocarcinoma	I	F	55
29	D2	Uterus	Endometrial adenocarcinoma	I	F	50
30	D3	Uterus	Endometrial adenocarcinoma	I	F	80
31	D4	Uterus	Adenosquamous carcinoma	-	F	65
32	D5	Uterus	Endometrial adenocarcinoma	II	F	62
33	D6	Uterus	Endometrial adenocarcinoma	II	F	58
34	D7	Uterus	Endometrial adenocarcinoma	I	F	57
35	D8	Uterus	Endometrial adenocarcinoma	II	F	51



PROTEIN

BIOTECHNOLOGIES

36	D9	Uterus	Endometrial adenocarcinoma	I	F	57
37	E1	Uterus	Endometrial adenocarcinoma	II	F	67
38	E2	Uterus	Endometrial adenocarcinoma	II	F	45
39	E3	Uterus	Endometrial adenocarcinoma	III	F	62
40	E4	Uterus	Endometrial adenocarcinoma	II	F	52
41	E5	Uterus	Endometrial adenocarcinoma	II	F	55
42	E6	Uterus	Endometrial adenocarcinoma	II	F	72
43	E7	Uterus	Endometrial adenocarcinoma	II	F	62
44	E8	Uterus	Adenosquamous carcinoma	-	F	49
45	E9	Uterus	Endometrial adenocarcinoma	II	F	51
46	F1	Uterus	Endometrial adenocarcinoma	I	F	62
47	F2	Uterus	Endometrial adenocarcinoma	III	F	53
48	F3	Uterus	Endometrial adenocarcinoma	III	F	48
49	F4	Uterus	Endometrial adenocarcinoma	II	F	39
50	F5	Uterus	Endometrial papillary adenocarcinoma	I	F	60
51	F6	Uterus	Endometrial adenocarcinoma	I	F	68
52	F7	Uterus	Endometrial adenocarcinoma	II	F	47
53	F8	Uterus	Endometrial adenocarcinoma	II	F	52
54	F9	Uterus	Endometrial adenocarcinoma	II	F	56
55	G1	Uterus	Endometrial adenocarcinoma	III	F	63
56	G2	Uterus	Endometrial adenocarcinoma	II	F	54
57	G3	Uterus	Endometrial adenocarcinoma	II	F	72
58	G4	Uterus	Endometrial adenocarcinoma	II	F	50
59	G5	Uterus	Endometrial adenocarcinoma	III	F	78
60	G6	Uterus	Endometrial adenocarcinoma	I	F	56
61	G7	Uterus	Secretory endometrium	-	F	48
62	G8	Uterus	Proliferative endometrium	-	F	56
63	G9	Uterus	Smooth muscle, blood vessel	-	F	61

For Research Use Only