



PROTEIN
BIOTECHNOLOGIES

HUMAN CERVIX CANCER TISSUE ARRAY

Catalog Number: TMA-071

Each dot represents a DISEASED tissue spot from one individual specimen that was selected and pathologically confirmed

Cases: 61 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** HPV confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Uterine cervix	Squamous cell carcinoma	II	F	62
02	A2	Uterine cervix	Inflammation tissue	-	F	61
03	A3	Uterine cervix	Squamous cell carcinoma	III	F	47
04	A4	Uterine cervix	Squamous cell carcinoma	II	F	51
05	A5	Uterine cervix	Squamous cell carcinoma	II	F	56
06	A6	Uterine cervix	Squamous cell carcinoma	II	F	30
07	A7	Uterine cervix	Squamous cell carcinoma	III	F	30
08	A8	Uterine cervix	Squamous cell carcinoma	III	F	39
09	A9	Uterine cervix	Smooth muscle	-	F	30
10	B1	Uterine cervix	Squamous cell carcinoma	II	F	45
11	B2	Uterine cervix	Squamous cell carcinoma	II	F	39
12	B3	Uterine cervix	Squamous cell carcinoma	II	F	35
13	B4	Uterine cervix	Squamous cell carcinoma	III	F	48
14	B5	Uterine cervix	Squamous cell carcinoma	II	F	45
15	B6	Uterine cervix	Squamous cell carcinoma	III	F	48
16	B7	Uterine cervix	Squamous cell carcinoma	II	F	43
17	B8	Uterine cervix	Squamous cell carcinoma	III	F	64
18	B9	Uterine cervix	Squamous cell carcinoma	II	F	37
19	C1	Uterine cervix	Squamous cell carcinoma	II	F	44
20	C2	Uterine cervix	Squamous cell carcinoma	II	F	49
21	C3	Uterine cervix	Squamous cell carcinoma	II	F	63
22	C4	Uterine cervix	Squamous cell carcinoma	III	F	42
23	C5	Uterine cervix	Squamous cell carcinoma	III	F	62
24	C6	Uterine cervix	Squamous cell carcinoma	II	F	27
25	C7	Uterine cervix	Mucosae tissue	-	F	39
26	C8	Uterine cervix	Squamous cell carcinoma	II	F	37
27	C9	Uterine cervix	Squamous cell carcinoma	II	F	43
28	D1	Uterine cervix	Squamous cell carcinoma	II	F	38
29	D2	Uterine cervix	Squamous cell carcinoma	III	F	42
30	D3	Uterine cervix	Squamous cell carcinoma	III	F	34
31	D4	Uterine cervix	Squamous cell carcinoma	II	F	56
32	D5	Uterine cervix	Squamous cell carcinoma	II	F	48
33	D6	Uterine cervix	Squamous cell carcinoma	II	F	54
34	D7	Uterine cervix	Squamous cell carcinoma	II	F	44
35	D8	Uterine cervix	Squamous cell carcinoma	II	F	38



PROTEIN

BIOTECHNOLOGIES

36	D9	Uterine cervix	Squamous cell carcinoma	III	F	55
37	E1	Uterine cervix	Squamous cell carcinoma	II	F	50
38	E2	Uterine cervix	Squamous cell carcinoma	II	F	54
39	E3	Uterine cervix	Squamous cell carcinoma	II	F	44
40	E4	Uterine cervix	Squamous cell carcinoma	II	F	41
41	E5	Uterine cervix	Squamous cell carcinoma	III	F	50
42	E6	Uterine cervix	Squamous cell carcinoma	III	F	61
43	E7	Uterine cervix	Squamous cell carcinoma	II	F	54
44	E8	Uterine cervix	Squamous cell carcinoma	II	F	43
45	E9	Uterine cervix	Squamous cell carcinoma	II	F	57
46	F1	Uterine cervix	Squamous cell carcinoma	II	F	32
47	F2	Uterine cervix	Squamous cell carcinoma	III	F	52
48	F3	Uterine cervix	Squamous cell carcinoma	III	F	38
49	F4	Uterine cervix	Squamous cell carcinoma	II	F	30
50	F5	Uterine cervix	Squamous cell carcinoma	II	F	44
51	F6	Uterine cervix	Squamous cell carcinoma	III	F	52
52	F7	Uterine cervix	Squamous cell carcinoma	III	F	38
53	F8	Uterine cervix	Adenosquamous carcinoma	II	F	36
54	F9	Uterine cervix	Adenocarcinoma	I	F	38
55	G1	Uterine cervix	Squamous cell carcinoma	II	F	55
56	G2	Uterine cervix	Squamous cell carcinoma	II	F	45
57	G3	Uterine cervix	Squamous cell carcinoma	II	F	33
58	G4	Uterine cervix	Squamous cell carcinoma	II	F	51
59	G5	Uterine cervix	Squamous cell carcinoma	III	F	56
60	G6	Uterine cervix	Squamous cell carcinoma	III	F	56
61	G7	Uterine cervix	Normal cervical tissue	-	F	39
62	G8	Uterine cervix	Normal cervical tissue	-	F	56
63	G9	Uterine cervix	Normal cervical tissue	-	F	55

For Research Use Only