



PROTEIN
BIOTECHNOLOGIES

HUMAN THYROID CANCER TISSUE ARRAY

Catalog Number: TMA-078

Each dot represents a DISEASED tissue spot from one individual specimen that was selected and pathologically confirmed

Cases: 80 **Cores:** 80 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** Thyroglobulin(TG) antibody confirmed

Layout: 8 x 10

	1	2	3	4	5	6	7	8	9	10
A	1	2	3	4	5	6	7	8	9	10
B	11	12	13	14	15	16	17	18	19	20
C	21	22	23	24	25	26	27	28	29	30
D	31	32	33	34	35	36	37	38	39	40
E	41	42	43	44	45	46	47	48	49	50
F	51	52	53	54	55	56	57	58	59	60
G	61	62	63	64	65	66	67	68	69	70
H	71	72	73	74	75	76	77	78	79	80

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Thyroid	Papillary carcinoma	-	F	56
02	A2	Thyroid	Papillary carcinoma	-	M	38
03	A3	Thyroid	Papillary carcinoma	-	F	30
04	A4	Thyroid	Papillary carcinoma	-	F	38
05	A5	Thyroid	Papillary carcinoma	-	F	38
06	A6	Thyroid	Lymphocytic cell infiltrating thyroid tissue	-	F	42
07	A7	Thyroid	Follicular papillary carcinoma	-	F	40
08	A8	Thyroid	Papillary carcinoma	-	F	18
09	A9	Thyroid	Papillary carcinoma	-	F	64
10	A10	Thyroid	Papillary carcinoma	-	M	33
11	B1	Thyroid	Papillary carcinoma	-	F	53
12	B2	Thyroid	Papillary carcinoma	-	F	56
13	B3	Thyroid	Papillary carcinoma	-	M	52
14	B4	Thyroid	Papillary carcinoma	-	F	34
15	B5	Thyroid	Lymphocytic cell infiltrating thyroid tissue	-	F	33
16	B6	Thyroid	Papillary carcinoma	-	F	40
17	B7	Thyroid	Follicular carcinoma	-	F	58
18	B8	Thyroid	Papillary carcinoma	-	F	37
19	B9	Thyroid	Papillary carcinoma	-	M	39
20	B10	Thyroid	Papillary carcinoma	-	F	46
21	C1	Thyroid	Papillary carcinoma	-	F	25
22	C2	Thyroid	Papillary carcinoma	-	F	21
23	C3	Thyroid	Follicular carcinoma	-	F	36
24	C4	Thyroid	Follicular carcinoma	-	M	56
25	C5	Thyroid	Follicular carcinoma	-	F	55
26	C6	Thyroid	Papillary carcinoma	-	F	68
27	C7	Thyroid	Follicular carcinoma	III	F	49
28	C8	Thyroid	Follicular carcinoma	-	F	70
29	C9	Thyroid	Follicular carcinoma	-	M	34
30	C10	Thyroid	Papillary carcinoma	-	F	79
31	D1	Thyroid	Papillary carcinoma	-	F	60
32	D2	Thyroid	Follicular carcinoma	-	F	60
33	D3	Thyroid	Medullary carcinoma	-	F	72



PROTEIN

BIOTECHNOLOGIES

34	D4	Thyroid	Medullary carcinoma	-	F	63
35	D5	Thyroid	Medullary carcinoma	-	F	71
36	D6	Thyroid	Medullary carcinoma	-	M	63
37	D7	Thyroid	Medullary carcinoma	-	M	51
38	D8	Thyroid	Medullary carcinoma	-	F	40
39	D9	Thyroid	Tall cell papillary carcinoma	-	F	52
40	D10	Thyroid	Tall cell papillary carcinoma	-	F	60
41	E1	Thyroid	Follicular adenoma	-	M	62
42	E2	Thyroid	Follicular adenoma	-	M	26
43	E3	Thyroid	Clear cell type adenoma	-	F	50
44	E4	Thyroid	Follicular adenoma	-	F	50
45	E5	Thyroid	Follicular adenoma	-	F	52
46	E6	Thyroid	Follicular adenoma	-	F	35
47	E7	Thyroid	Follicular adenoma	-	F	31
48	E8	Thyroid	Follicular adenoma	-	F	21
49	E9	Thyroid	Follicular adenoma	-	F	56
50	E10	Thyroid	Follicular adenoma	-	F	31
51	F1	Thyroid	Follicular adenoma	-	F	31
52	F2	Thyroid	Colloid adenoma	-	F	36
53	F3	Thyroid	Colloid adenoma	-	F	21
54	F4	Thyroid	Colloid adenoma	-	F	42
55	F5	Thyroid	Colloid adenoma	-	F	32
56	F6	Thyroid	Colloid adenoma	-	F	35
57	F7	Thyroid	Colloid adenoma	-	F	48
58	F8	Thyroid	Colloid adenoma	-	F	48
59	F9	Thyroid	Colloid adenoma	-	M	35
60	F10	Thyroid	Colloid adenoma	-	F	52
61	G1	Thyroid	Colloid adenoma	-	F	26
62	G2	Thyroid	Coagulated blood	-	F	45
63	G3	Thyroid	Colloid adenoma	-	F	38
64	G4	Thyroid	Colloid adenoma	-	F	46
65	G5	Thyroid	Colloid adenoma	-	F	70
66	G6	Thyroid	Colloid adenoma	-	M	55
67	G7	Thyroid	Embryonal adenoma	-	F	46
68	G8	Thyroid	Embryonal adenoma	-	F	25
69	G9	Thyroid	Embryonal adenoma	-	M	46
70	G10	Thyroid	Embryonal adenoma	-	F	39
71	H1	Thyroid	Embryonal adenoma	-	F	38
72	H2	Thyroid	Embryonal adenoma	-	F	20
73	H3	Thyroid	Embryonal adenoma	-	F	46
74	H4	Thyroid	Embryonal adenoma	-	F	52
75	H5	Thyroid	Embryonal adenoma	-	F	19
76	H6	Thyroid	Embryonal adenoma	-	M	42
77	H7	Thyroid	Oxyphilic adenoma	-	F	56
78	H8	Thyroid	Oxyphilic adenoma	-	F	37
79	H9	Thyroid	Oxyphilic adenoma	-	F	40
80	H10	Thyroid	Atypical adenoma	-	F	33

For Research Use Only