



PROTEIN
BIOTECHNOLOGIES

HUMAN STOMACH CANCER TISSUE ARRAY

Catalog Number: TMA-094

Each two (2) dots from two different tissue spots represents one single specimen that was pathologically selected and confirmed (Gastric ulcer tissue array)

Cases: 14 **Cores:** 24 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** CK confirmed

Layout: 4 x 6

	1	2	3	4	5	6
A	1	2	3	4	5	6
B	7	8	9	10	11	12
C	13	14	15	16	17	18
D	19	20	21	22	23	24

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Stomach	Chronic gastric ulcer	-	M	27
02	A2	Stomach	Chronic gastric ulcer	-	M	27
03	A3	Stomach	Chronic gastric ulcer	-	M	66
04	A4	Stomach	Chronic gastric ulcer	-	M	66
05	A5	Stomach	Chronic gastric ulcer	-	M	61
06	A6	Stomach	Chronic gastric ulcer	-	M	61
07	B1	Stomach	Fibrous tissue blood vessel and inflammation"	-	F	52
08	B2	Stomach	Chronic gastric ulcer	-	F	52
09	B3	Stomach	Chronic gastric ulcer	-	M	61
10	B4	Stomach	Smooth muscle blood vessel and fibrous tissue"	-	M	61
11	B5	Stomach	Chronic gastric ulcer	-	M	57
12	B6	Stomach	Chronic gastric ulcer	-	M	57
13	C1	Stomach	Chronic gastric ulcer	-	M	36
14	C2	Stomach	Chronic gastric ulcer	-	M	36
15	C3	Stomach	Chronic gastric ulcer	-	M	65
16	C4	Stomach	Chronic gastric ulcer	-	M	65
17	C5	Stomach	Chronic gastric ulcer	-	M	80
18	C6	Stomach	Chronic gastric ulcer	-	M	80
19	D1	Stomach	Chronic gastric ulcer	-	M	56
20	D2	Stomach	Chronic gastric ulcer	-	M	56
21	D3	Stomach	Mucosa	-	M	56
22	D4	Stomach	Mucosa	-	M	62
23	D5	Stomach	Mucosa	-	M	42
24	D6	Stomach	Mucosa	-	M	56

For Research Use Only