



PROTEIN
BIOTECHNOLOGIES

HUMAN STOMACH CANCER TISSUE ARRAY

Catalog Number: TMA-095

Each dot represents a DISEASED tissue spot from one individual specimen that was selected and pathologically confirmed (Stomach Carcinoma Tissue Array)

Cases: 0 **Cores:** 72 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** CgA confirmed

Layout: 8 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63
H	64	65	66	67	68	69	70	71	72

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Stomach	Adenocarcinoma	III	M	54
02	A2	Stomach	Adenocarcinoma	III	M	37
03	A3	Stomach	Chronic inflammation of mucosa	-	F	57
04	A4	Stomach	Adenocarcinoma	III	M	66
05	A5	Stomach	Adenocarcinoma	II	F	41
06	A6	Stomach	Adenocarcinoma	III	F	65
07	A7	Stomach	Adenocarcinoma	III	F	54
08	A8	Stomach	Adenocarcinoma	II	M	69
09	A9	Stomach	Adenocarcinoma	III	M	72
10	B1	Stomach	Adenocarcinoma	III	M	69
11	B2	Stomach	Adenocarcinoma	III	M	63
12	B3	Stomach	Squamous cell carcinoma	I	M	50
13	B4	Stomach	Carcinoid	-	M	62
14	B5	Stomach	Adenocarcinoma	III	M	56
15	B6	Stomach	Adenocarcinoma	II	M	75
16	B7	Stomach	Adenocarcinoma	III	M	81
17	B8	Stomach	Adenocarcinoma	II	M	66
18	B9	Stomach	Adenocarcinoma	II	M	64
19	C1	Stomach	Mucinous adenocarcinoma	III	M	55
20	C2	Stomach	Squamous cell carcinoma	II	M	65
21	C3	Stomach	Carcinoid	-	M	71
22	C4	Stomach	Adenosquamous carcinoma	-	F	68
23	C5	Stomach	Adenocarcinoma	III	M	53
24	C6	Stomach	Carcinoid	-	M	61
25	C7	Stomach	Squamous cell carcinoma	III	F	48
26	C8	Stomach	Chronic inflammation of mucosa	-	M	58
27	C9	Stomach	Adenocarcinoma	III	M	47
28	D1	Stomach	Adenocarcinoma	III	M	62
29	D2	Stomach	Adenocarcinoma	III	M	72
30	D3	Stomach	Adenocarcinoma	II	F	60
31	D4	Stomach	Adenocarcinoma	II	M	65
32	D5	Stomach	Adenocarcinoma	II	M	48
33	D6	Stomach	Adenocarcinoma	I	M	56
34	D7	Stomach	Adenocarcinoma	III	M	40



PROTEIN

BIOTECHNOLOGIES

35	D8	Stomach	Adenocarcinoma	II	M	65
36	D9	Stomach	Adenocarcinoma	III	M	68
37	E1	Stomach	Adenocarcinoma	III	F	35
38	E2	Stomach	Adenocarcinoma	III	M	63
39	E3	Stomach	Adenocarcinoma	III	F	31
40	E4	Stomach	Adenocarcinoma	II	M	60
41	E5	Stomach	Adenocarcinoma	III	M	66
42	E6	Stomach	Adenocarcinoma	III	M	60
43	E7	Stomach	Smooth muscle	-	F	31
44	E8	Stomach	Adenocarcinoma	I	M	61
45	E9	Stomach	Adenocarcinoma	II	M	56
46	F1	Stomach	Adenocarcinoma	III	F	28
47	F2	Stomach	Adenocarcinoma	III	M	52
48	F3	Stomach	Adenocarcinoma	II	F	61
49	F4	Stomach	Adenocarcinoma	II	M	58
50	F5	Stomach	Adenocarcinoma	III	M	77
51	F6	Stomach	Adenocarcinoma	III	F	80
52	F7	Stomach	Adenocarcinoma	II	M	72
53	F8	Stomach	Adenocarcinoma	II	M	44
54	F9	Stomach	Mucinous adenocarcinoma	III	F	52
55	G1	Stomach	Adenocarcinoma	II	M	56
56	G2	Stomach	Adenocarcinoma	II	M	44
57	G3	Stomach	Adenocarcinoma	II	M	35
58	G4	Stomach	Adenocarcinoma	III	M	48
59	G5	Stomach	Adenocarcinoma	II	M	47
60	G6	Stomach	Adenocarcinoma	II	M	75
61	G7	Stomach	Adenocarcinoma	III	M	61
62	G8	Stomach	Adenocarcinoma	II	M	42
63	G9	Stomach	Adenocarcinoma	II	M	65
64	H1	Stomach	Adenocarcinoma	III	M	66
65	H2	Stomach	Adenocarcinoma	III	M	60
66	H3	Stomach	Adenocarcinoma	II	F	50
67	H4	Stomach	Adenocarcinoma	III	M	65
68	H5	Stomach	Adenocarcinoma	III	M	64
69	H6	Stomach	Mucinous adenocarcinoma	III	M	72
70	H7	Stomach	Mucosa	-	M	71
71	H8	Stomach	Mucosa	-	M	54
72	H9	Stomach	Mucosa	-	M	53

For Research Use Only