



PROTEIN
BIOTECHNOLOGIES

HUMAN STOMACH CANCER TISSUE ARRAY

Catalog Number: TMA-097

This is a COMBINATION ARRAY PANEL from diseased, normal and adjacent (1.5cm) to specimen that was selected and pathologically confirmed (stomach adenocarcinoma tissue array)

Cases: 21 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** L-CK confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Stomach	Mucosae chronic inflammation	-	F	55
02	A2	Stomach	Mucosae chronic inflammation	-	F	55
03	A3	Stomach	Mucosae chronic inflammation	-	F	55
04	A4	Stomach	Fatty tissue and blood vessel	-	M	64
05	A5	Stomach	Adenocarcinoma	II	M	64
06	A6	Stomach	Mucosae chronic inflammation	-	M	64
07	A7	Stomach	Mucosae	-	M	43
08	A8	Stomach	Adenocarcinoma	III	M	43
09	A9	Stomach	Mucosae chronic inflammation	-	M	43
10	B1	Stomach	Mucosae	-	M	71
11	B2	Stomach	Mucinous adenocarcinoma	III	M	71
12	B3	Stomach	Mucosae chronic inflammation	-	M	71
13	B4	Stomach	Mucosae	-	M	53
14	B5	Stomach	Smooth muscle and blood vessel	-	M	53
15	B6	Stomach	Mucosae	-	M	53
16	B7	Stomach	Mucosae	-	M	30
17	B8	Stomach	Adenocarcinoma	III	M	30
18	B9	Stomach	Mucosae chronic inflammation	-	M	30
19	C1	Stomach	Mucosae	-	M	69
20	C2	Stomach	Fatty tissue and blood vessel	-	M	69
21	C3	Stomach	Adenocarcinoma	III	M	69
22	C4	Stomach	Mucosae	-	M	33
23	C5	Stomach	Mucinous adenocarcinoma	II	M	33
24	C6	Stomach	A few mucosae chronic inflammation	-	M	33
25	C7	Stomach	Smooth muscle and fatty tissue	-	F	60
26	C8	Stomach	Adenocarcinoma	III	F	60
27	C9	Stomach	Mucosae chronic inflammation	-	F	60
28	D1	Stomach	Mucosae chronic inflammation	-	M	58
29	D2	Stomach	Mucosae chronic inflammation	-	M	58
30	D3	Stomach	Mucosae chronic inflammation	-	M	58
31	D4	Stomach	Smooth muscle	-	M	47
32	D5	Stomach	Adenocarcinoma	III	M	47
33	D6	Stomach	Fibrous tissue and blood vessel	-	M	47
34	D7	Stomach	Fibrofatty tissue	-	M	65
35	D8	Stomach	Adenocarcinoma	III	M	65



PROTEIN

BIOTECHNOLOGIES

36	D9	Stomach	Fatty tissue and blood vessel	-	M	65
37	E1	Stomach	Smooth muscle	-	M	67
38	E2	Stomach	Adenocarcinoma	III	M	67
39	E3	Stomach	Chronic inflammation of esophagus mucosae	-	M	67
40	E4	Stomach	Mucosae	-	F	73
41	E5	Stomach	Mucinous adenocarcinoma	III	F	73
42	E6	Stomach	Mucosae chronic inflammation	-	F	73
43	E7	Stomach	Mucosae	-	M	48
44	E8	Stomach	Mucinous adenocarcinoma	I	M	48
45	E9	Stomach	Mucosae chronic inflammation	-	M	48
46	F1	Stomach	Mucosae	-	M	60
47	F2	Stomach	Adenocarcinoma	I	M	60
48	F3	Stomach	Mucosae chronic inflammation	-	M	60
49	F4	Stomach	Mucosae	-	M	45
50	F5	Stomach	Adenocarcinoma	I	M	45
51	F6	Stomach	Mild atypical hyperplasia	-	M	45
52	F7	Stomach	Mucosae chronic inflammation	-	M	56
53	F8	Stomach	Adenocarcinoma	III	M	56
54	F9	Stomach	Mucosae chronic inflammation	-	M	56
55	G1	Stomach	Smooth muscle and blood vessel	-	F	31
56	G2	Stomach	Adenocarcinoma	III	F	31
57	G3	Stomach	Mucosae chronic inflammation	-	F	31
58	G4	Stomach	Mucosae	-	F	43
59	G5	Stomach	Mucinous adenocarcinoma	III	F	43
60	G6	Stomach	Mucosae chronic inflammation	-	F	43
61	G7	Stomach	Mucosae	-	F	49
62	G8	Stomach	Adenocarcinoma	II	F	49
63	G9	Stomach	Fibrous and blood vessel	-	F	49

For Research Use Only