



**PROTEIN**  
BIOTECHNOLOGIES

**HUMAN STOMACH CANCER TISSUE ARRAY**

**Catalog Number:** TMA-098

This is a COMBINATION ARRAY PANEL from diseased, normal and adjacent (1.5cm) to specimen that was selected and pathologically confirmed (stomach adenocarcinoma tissue array)

**Cases:** 21    **Cores:** 63    **Diameter:** 1.5mm    **Thickness:** 5 µm    **Standard IHC:** CK confirmed

**Layout:** 7 x 9

	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>
<b>A</b>	1	2	3	4	5	6	7	8	9
<b>B</b>	10	11	12	13	14	15	16	17	18
<b>C</b>	19	20	21	22	23	24	25	26	27
<b>D</b>	28	29	30	31	32	33	34	35	36
<b>E</b>	37	38	39	40	41	42	43	44	45
<b>F</b>	46	47	48	49	50	51	52	53	54
<b>G</b>	55	56	57	58	59	60	61	62	63

<b>Core No.</b>	<b>Location</b>	<b>Tissue</b>	<b>Diagnosis</b>	<b>Grade</b>	<b>Sex</b>	<b>Age</b>
01	A1	Stomach	Adenocarcinoma	III	M	58
02	A2	Stomach	Loose fibrous tissue and blood vessel	-	M	58
03	A3	Stomach	Mild chronic inflammation of mucosa	-	M	58
04	A4	Stomach	Adenocarcinoma, mainly is mucous cell carcinoma	III	M	55
05	A5	Stomach	Moderate atypical hyperplasia of gland	-	M	55
06	A6	Stomach	A few glands of mucosa	-	M	55
07	A7	Stomach	Adenocarcinoma	III	M	53
08	A8	Stomach	Chronic inflammation of mucosa	-	M	53
09	A9	Stomach	Chronic inflammation of mucosa	-	M	53
10	B1	Stomach	Adenocarcinoma	III	M	59
11	B2	Stomach	Loose fibrous tissue	-	M	59
12	B3	Stomach	Chronic inflammation of mucosa	-	M	59
13	B4	Stomach	Adenocarcinoma	III	F	50
14	B5	Stomach	Mild hyperplasia of gland	-	F	50
15	B6	Stomach	Mucosa tissue	-	F	50
16	B7	Stomach	Adenocarcinoma	III	M	44
17	B8	Stomach	Loose fibrous tissue	-	M	44
18	B9	Stomach	Mild chronic inflammation of mucosa	-	M	44
19	C1	Stomach	Adenocarcinoma	II	M	64
20	C2	Stomach	Moderate atypical hyperplasia of gland	-	M	64
21	C3	Stomach	Mild chronic inflammation of mucosa	-	M	64
22	C4	Stomach	Adenocarcinoma	II	M	69
23	C5	Stomach	Moderate atypical hyperplasia of gland	-	M	69
24	C6	Stomach	Mucosa tissue	-	M	69
25	C7	Stomach	Adenocarcinoma	I	F	59
26	C8	Stomach	Mild atypical hyperplasia of gland	-	F	59
27	C9	Stomach	Normal mucosa tissue	-	F	59
28	D1	Stomach	Adenocarcinoma	III	M	63
29	D2	Stomach	Mild hyperplasia of gland of a small area of mucosa	-	M	63
30	D3	Stomach	Superficial gastritis	-	M	63
31	D4	Stomach	Adenocarcinoma	III	M	62
32	D5	Stomach	Superficial gastritis	-	M	62
33	D6	Stomach	Glands of mucosa	-	M	62
34	D7	Stomach	Adenocarcinoma	III	M	55
35	D8	Stomach	Mild hyperplasia of a few glands	-	M	55



**PROTEIN**  
BIOTECHNOLOGIES

36	D9	Stomach	Mucosa tissue	-	M	55
37	E1	Stomach	Adenocarcinoma	III	M	45
38	E2	Stomach	Chronic atrophic gastritis	-	M	45
39	E3	Stomach	Chronic inflammation of mucosa	-	M	45
40	E4	Stomach	Adenocarcinoma	III	M	41
41	E5	Stomach	Mild atypical hyperplasia of gland	-	M	41
42	E6	Stomach	Chronic superficial gastritis	-	M	41
43	E7	Stomach	Adenocarcinoma	-	F	43
44	E8	Stomach	Smooth muscle and fibrous tissue	-	F	43
45	E9	Stomach	A few glands of mucosa	-	F	43
46	F1	Stomach	Adenocarcinoma	II	M	47
47	F2	Stomach	Chronic inflammation of mucosa, partly mild atypical hyperplasia of gland	-	M	47
48	F3	Stomach	Chronic inflammation of mucosa	-	M	47
49	F4	Stomach	Adenocarcinoma	III	M	51
50	F5	Stomach	Mild chronic inflammation of mucosa	-	M	51
51	F6	Stomach	Mild chronic inflammation of mucosa	-	M	51
52	F7	Stomach	Adenocarcinoma	III	M	67
53	F8	Stomach	Adenocarcinoma	I	M	67
54	F9	Stomach	Glands of mucosa	-	M	67
55	G1	Stomach	Adenocarcinoma	III	M	55
56	G2	Stomach	Chronic inflammation of mucosa	-	M	55
57	G3	Stomach	Chronic superficial gastritis	-	M	55
58	G4	Stomach	Adenocarcinoma	III	M	67
59	G5	Stomach	Mild chronic inflammation of mucosa	-	M	67
60	G6	Stomach	Mild chronic inflammation of mucosa	-	M	67
61	G7	Stomach	Adenocarcinoma	I	M	71
62	G8	Stomach	Chronic inflammation of mucosa	-	M	71
63	G9	Stomach	Chronic inflammation of mucosa	-	M	71

**For Research Use Only**