



PROTEIN
BIOTECHNOLOGIES

HUMAN STOMACH CANCER TISSUE ARRAY

Catalog Number: TMA-102

Each dot represents a NORMAL tissue spot from one individual specimen that was selected and pathologically confirmed

Cases: 70 **Cores:** 72 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** Actin confirmed

Layout: 8 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63
H	64	65	66	67	68	69	70	71	72

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Stomach	Gastric mucosa	-	F	55
02	A2	Stomach	Gastric mucosa	-	M	59
03	A3	Stomach	Gastric mucosa	-	M	66
04	A4	Stomach	Gastric mucosa	-	M	54
05	A5	Stomach	Gastric mucosa	-	M	43
06	A6	Stomach	Gastric mucosa	-	F	65
07	A7	Stomach	Gastric mucosa	-	M	65
08	A8	Stomach	Smooth muscle and fatty tissue	-	F	54
09	A9	Stomach	Small area of mucosae	-	M	69
10	B1	Stomach	Gastric mucosa	-	F	65
11	B2	Stomach	Gastric mucosa	-	M	72
12	B3	Stomach	Gastric mucosa	-	M	44
13	B4	Stomach	Gastric mucosa	-	M	60
14	B5	Stomach	Gastric mucosa	-	M	60
15	B6	Stomach	Smooth muscle and fatty tissue	-	M	55
16	B7	Stomach	Gastric mucosa	-	M	56
17	B8	Stomach	Gastric mucosa	-	F	34
18	B9	Stomach	Smooth muscle	-	M	66
19	C1	Stomach	Smooth muscle and blood vessel	-	M	55
20	C2	Stomach	Gastric mucosa	-	M	68
21	C3	Stomach	Gastric mucosa	-	M	55
22	C4	Stomach	Blood vessel and fatty tissue	-	M	77
23	C5	Stomach	Gastric mucosa	-	M	71
24	C6	Stomach	Blood vessel and fatty tissue	-	M	53
25	C7	Stomach	Gastric mucosa	-	M	69
26	C8	Stomach	Gastric mucosa	-	M	33
27	C9	Stomach	Fibrofatty tissue	-	M	58
28	D1	Stomach	Fibrofatty tissue	-	M	47
29	D2	Stomach	Gastric mucosa	-	M	30
30	D3	Stomach	Adenocarcinoma infiltration	II	M	58
31	D4	Stomach	Fibrofatty tissue	-	F	60
32	D5	Stomach	Gastric mucosa	-	F	52
33	D6	Stomach	Gastric mucosa	-	M	48
34	D7	Stomach	Gastric mucosa	-	M	56



PROTEIN

BIOTECHNOLOGIES

35	D8	Stomach	Blood vessel and fatty tissue	-	M	57
36	D9	Stomach	Gastric mucosa	-	M	65
37	E1	Stomach	Gastric mucosa	-	F	35
38	E2	Stomach	Gastric mucosa	-	M	68
39	E3	Stomach	Gastric mucosa	-	M	60
40	E4	Stomach	Gastric mucosa	-	F	55
41	E5	Stomach	Small area of mucosae	-	M	45
42	E6	Stomach	Gastric mucosa	-	M	52
43	E7	Stomach	Gastric mucosa	-	F	57
44	E8	Stomach	Gastric mucosa	-	F	50
45	E9	Stomach	Smooth muscle and blood vessel	-	F	31
46	F1	Stomach	Gastric mucosa	-	F	43
47	F2	Stomach	Blood vessel and fatty tissue	-	M	43
48	F3	Stomach	Gastric mucosa	-	M	60
49	F4	Stomach	Gastric mucosa	-	M	50
50	F5	Stomach	Gastric mucosa	-	M	61
51	F6	Stomach	Gastric mucosa	-	M	56
52	F7	Stomach	Gastric mucosa	-	M	41
53	F8	Stomach	Gastric mucosa	-	M	66
54	F9	Stomach	Gastric mucosa	-	F	61
55	G1	Stomach	Gastric mucosa	-	M	58
56	G2	Stomach	Gastric mucosa	-	F	63
57	G3	Stomach	Gastric mucosa	-	F	80
58	G4	Stomach	Small area of mucosae	-	M	72
59	G5	Stomach	Gastric mucosa	-	M	44
60	G6	Stomach	Gastric mucosa	-	M	62
61	G7	Stomach	Gastric mucosa	-	M	44
62	G8	Stomach	Gastric mucosa	-	M	47
63	G9	Stomach	Gastric mucosa	-	M	67
64	H1	Stomach	Gastric mucosa	-	M	55
65	H2	Stomach	Gastric mucosa	-	F	46
66	H3	Stomach	Gastric mucosa	-	M	63
67	H4	Stomach	Gastric mucosa	-	M	59
68	H5	Stomach	Gastric mucosa	-	M	64
69	H6	Stomach	Gastric mucosa	-	M	54
70	H7	Stomach	Adenocarcinoma	III	F	55
71	H8	Stomach	Smooth muscle	-	M	57
72	H9	Stomach	Carcinoma cells infiltration	-	M	73

For Research Use Only