



# PROTEIN BIOTECHNOLOGIES

## HUMAN STOMACH CANCER TISSUE ARRAY

**Catalog Number:** TMA-105

Each three (3) dots from three different tissue spots represents one single specimen selected and pathologically confirmed. (Stomach Adenocarcinoma tissue array)

**Cases:** 22    **Cores:** 63    **Diameter:** 1.5mm    **Thickness:** 5 µm    **Standard IHC:** L-CK confirmed

**Layout:** 7 x 9

|          | <b>1</b> | <b>2</b> | <b>3</b> | <b>4</b> | <b>5</b> | <b>6</b> | <b>7</b> | <b>8</b> | <b>9</b> |
|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|
| <b>A</b> | 1        | 2        | 3        | 4        | 5        | 6        | 7        | 8        | 9        |
| <b>B</b> | 10       | 11       | 12       | 13       | 14       | 15       | 16       | 17       | 18       |
| <b>C</b> | 19       | 20       | 21       | 22       | 23       | 24       | 25       | 26       | 27       |
| <b>D</b> | 28       | 29       | 30       | 31       | 32       | 33       | 34       | 35       | 36       |
| <b>E</b> | 37       | 38       | 39       | 40       | 41       | 42       | 43       | 44       | 45       |
| <b>F</b> | 46       | 47       | 48       | 49       | 50       | 51       | 52       | 53       | 54       |
| <b>G</b> | 55       | 56       | 57       | 58       | 59       | 60       | 61       | 62       | 63       |

| <b>Core No.</b> | <b>Location</b> | <b>Tissue</b> | <b>Diagnosis</b>    | <b>Grade</b> | <b>Sex</b> | <b>Age</b> |
|-----------------|-----------------|---------------|---------------------|--------------|------------|------------|
| 01              | A1              | Stomach       | Adenocarcinoma      | III          | F          | 55         |
| 02              | A2              | Stomach       | Adenocarcinoma      | III          | F          | 55         |
| 03              | A3              | Stomach       | Adenocarcinoma      | III          | F          | 55         |
| 04              | A4              | Stomach       | Inflammatory tissue | -            | M          | 57         |
| 05              | A5              | Stomach       | Inflammatory tissue | -            | M          | 57         |
| 06              | A6              | Stomach       | Inflammatory tissue | -            | M          | 57         |
| 07              | A7              | Stomach       | Adenocarcinoma      | III          | M          | 55         |
| 08              | A8              | Stomach       | Adenocarcinoma      | III          | M          | 55         |
| 09              | A9              | Stomach       | Adenocarcinoma      | III          | M          | 55         |
| 10              | B1              | Stomach       | Adenocarcinoma      | III          | F          | 59         |
| 11              | B2              | Stomach       | Adenocarcinoma      | III          | F          | 59         |
| 12              | B3              | Stomach       | Adenocarcinoma      | III          | F          | 59         |
| 13              | B4              | Stomach       | Adenocarcinoma      | III          | M          | 69         |
| 14              | B5              | Stomach       | Adenocarcinoma      | III          | M          | 69         |
| 15              | B6              | Stomach       | Adenocarcinoma      | III          | M          | 69         |
| 16              | B7              | Stomach       | Adenocarcinoma      | III          | F          | 63         |
| 17              | B8              | Stomach       | Adenocarcinoma      | III          | F          | 63         |
| 18              | B9              | Stomach       | Adenocarcinoma      | III          | F          | 63         |
| 19              | C1              | Stomach       | Adenocarcinoma      | III          | M          | 64         |
| 20              | C2              | Stomach       | Adenocarcinoma      | III          | M          | 64         |
| 21              | C3              | Stomach       | Adenocarcinoma      | III          | M          | 64         |
| 22              | C4              | Stomach       | Adenocarcinoma      | III          | M          | 44         |
| 23              | C5              | Stomach       | Adenocarcinoma      | III          | M          | 44         |
| 24              | C6              | Stomach       | Adenocarcinoma      | III          | M          | 44         |
| 25              | C7              | Stomach       | Adenocarcinoma      | III          | M          | 56         |
| 26              | C8              | Stomach       | Adenocarcinoma      | III          | M          | 56         |
| 27              | C9              | Stomach       | Adenocarcinoma      | III          | M          | 56         |
| 28              | D1              | Stomach       | Adenocarcinoma      | III          | F          | 42         |
| 29              | D2              | Stomach       | Adenocarcinoma      | III          | F          | 42         |
| 30              | D3              | Stomach       | Adenocarcinoma      | III          | F          | 42         |
| 31              | D4              | Stomach       | Adenocarcinoma      | III          | M          | 55         |
| 32              | D5              | Stomach       | Adenocarcinoma      | III          | M          | 55         |
| 33              | D6              | Stomach       | Adenocarcinoma      | III          | M          | 55         |
| 34              | D7              | Stomach       | Adenocarcinoma      | III          | M          | 60         |
| 35              | D8              | Stomach       | Adenocarcinoma      | III          | M          | 60         |



# PROTEIN

BIOTECHNOLOGIES

|    |    |         |                |     |   |    |
|----|----|---------|----------------|-----|---|----|
| 36 | D9 | Stomach | Adenocarcinoma | III | M | 60 |
| 37 | E1 | Stomach | Adenocarcinoma | III | M | 50 |
| 38 | E2 | Stomach | Adenocarcinoma | III | M | 50 |
| 39 | E3 | Stomach | Adenocarcinoma | III | M | 50 |
| 40 | E4 | Stomach | Adenocarcinoma | III | M | 62 |
| 41 | E5 | Stomach | Adenocarcinoma | III | M | 62 |
| 42 | E6 | Stomach | Adenocarcinoma | III | M | 62 |
| 43 | E7 | Stomach | Adenocarcinoma | III | F | 65 |
| 44 | E8 | Stomach | Adenocarcinoma | III | F | 65 |
| 45 | E9 | Stomach | Adenocarcinoma | III | F | 65 |
| 46 | F1 | Stomach | Adenocarcinoma | III | M | 60 |
| 47 | F2 | Stomach | Adenocarcinoma | III | M | 60 |
| 48 | F3 | Stomach | Adenocarcinoma | III | M | 60 |
| 49 | F4 | Stomach | Adenocarcinoma | III | M | 69 |
| 50 | F5 | Stomach | Adenocarcinoma | III | M | 69 |
| 51 | F6 | Stomach | Adenocarcinoma | III | M | 69 |
| 52 | F7 | Stomach | Adenocarcinoma | III | M | 67 |
| 53 | F8 | Stomach | Adenocarcinoma | III | M | 67 |
| 54 | F9 | Stomach | Adenocarcinoma | III | M | 67 |
| 55 | G1 | Stomach | Adenocarcinoma | III | M | 74 |
| 56 | G2 | Stomach | Adenocarcinoma | III | M | 74 |
| 57 | G3 | Stomach | Adenocarcinoma | III | M | 74 |
| 58 | G4 | Stomach | Adenocarcinoma | III | M | 57 |
| 59 | G5 | Stomach | Adenocarcinoma | III | M | 57 |
| 60 | G6 | Stomach | Adenocarcinoma | III | M | 57 |
| 61 | G7 | Stomach | Normal tissue  | -   | F | 59 |
| 62 | G8 | Stomach | Normal tissue  | -   | M | 55 |
| 63 | G9 | Stomach | Normal tissue  | -   | M | 65 |

**For Research Use Only**