



PROTEIN
BIOTECHNOLOGIES

HUMAN STOMACH CANCER TISSUE ARRAY

Catalog Number: TMA-106

Each three (3) dots from three different tissue spots represents one single specimen selected and pathologically confirmed. (Stomach Adenocarcinoma tissue array)

Cases: 22 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** L-CK confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Stomach	Adenocarcinoma	III	M	64
02	A2	Stomach	Adenocarcinoma	III	M	64
03	A3	Stomach	Adenocarcinoma	III	M	64
04	A4	Stomach	Adenocarcinoma	III	M	57
05	A5	Stomach	Adenocarcinoma	III	M	57
06	A6	Stomach	Adenocarcinoma	III	M	57
07	A7	Stomach	Adenocarcinoma	III	M	55
08	A8	Stomach	Adenocarcinoma	III	M	55
09	A9	Stomach	Adenocarcinoma	III	M	55
10	B1	Stomach	Adenocarcinoma	III	M	68
11	B2	Stomach	Adenocarcinoma	III	M	68
12	B3	Stomach	Adenocarcinoma	III	M	68
13	B4	Stomach	Adenocarcinoma	III	M	69
14	B5	Stomach	Adenocarcinoma	III	M	69
15	B6	Stomach	Adenocarcinoma	III	M	69
16	B7	Stomach	Adenocarcinoma	III	M	54
17	B8	Stomach	Adenocarcinoma	III	M	54
18	B9	Stomach	Adenocarcinoma	III	M	54
19	C1	Stomach	Adenocarcinoma	III	M	64
20	C2	Stomach	Adenocarcinoma	III	M	64
21	C3	Stomach	Adenocarcinoma	III	M	64
22	C4	Stomach	Adenocarcinoma	III	M	44
23	C5	Stomach	Adenocarcinoma	III	M	44
24	C6	Stomach	Adenocarcinoma	III	M	44
25	C7	Stomach	Adenocarcinoma	III	M	56
26	C8	Stomach	Adenocarcinoma	III	M	56
27	C9	Stomach	Smooth muscle and fibrosum with chronic inflammation	-	M	56
28	D1	Stomach	Adenocarcinoma	III	M	59
29	D2	Stomach	Chronic inflammation	-	M	59
30	D3	Stomach	Adenocarcinoma	III	M	59
31	D4	Stomach	Adenocarcinoma	III	M	94
32	D5	Stomach	Adenocarcinoma	III	M	94
33	D6	Stomach	Adenocarcinoma	III	M	94
34	D7	Stomach	Adenocarcinoma	III	M	60



PROTEIN

BIOTECHNOLOGIES

35	D8	Stomach	Adenocarcinoma	III	M	60
36	D9	Stomach	Adenocarcinoma	III	M	60
37	E1	Stomach	Adenocarcinoma	III	M	70
38	E2	Stomach	Adenocarcinoma	III	M	70
39	E3	Stomach	Adenocarcinoma	III	M	70
40	E4	Stomach	Adenocarcinoma	III	M	60
41	E5	Stomach	Adenocarcinoma	III	M	60
42	E6	Stomach	Adenocarcinoma	III	M	60
43	E7	Stomach	Adenocarcinoma	III	M	65
44	E8	Stomach	Adenocarcinoma	III	M	65
45	E9	Stomach	Adenocarcinoma	III	M	65
46	F1	Stomach	Adenocarcinoma	III	M	60
47	F2	Stomach	Adenocarcinoma	III	M	60
48	F3	Stomach	Adenocarcinoma	III	M	60
49	F4	Stomach	Adenocarcinoma	III	M	35
50	F5	Stomach	Adenocarcinoma	III	M	35
51	F6	Stomach	Adenocarcinoma	III	M	35
52	F7	Stomach	Adenocarcinoma	III	M	49
53	F8	Stomach	Adenocarcinoma	III	M	49
54	F9	Stomach	Adenocarcinoma	III	M	49
55	G1	Stomach	Adenocarcinoma	III	M	74
56	G2	Stomach	Adenocarcinoma	III	M	74
57	G3	Stomach	Adenocarcinoma	III	M	74
58	G4	Stomach	Adenocarcinoma	III	F	45
59	G5	Stomach	Adenocarcinoma	III	F	45
60	G6	Stomach	Adenocarcinoma	III	F	45
61	G7	Stomach	Chronic inflammation	-	M	47
62	G8	Stomach	Chronic inflammation	-	M	55
63	G9	Stomach	Mucosae	-	M	57

For Research Use Only