



PROTEIN
BIOTECHNOLOGIES

HUMAN LIVER CANCER TISSUE ARRAY

Catalog Number: TMA-123

This is a COMBINATION ARRAY PANEL from diseased, normal and adjacent (1.5cm) to specimen that was selected and pathologically confirmed (Liver carcinoma tissue array)

Cases: 42 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** PCNA confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Liver	Viral hepatitis	-	F	66
02	A2	Liver	Adenocarcinoma	II	F	66
03	A3	Liver	Margin of carcinoma	-	M	43
04	A4	Liver	Viral hepatitis	-	F	37
05	A5	Liver	Hepatocellular carcinoma	I	M	55
06	A6	Liver	Margin of carcinoma	-	M	55
07	A7	Liver	Viral hepatitis	-	F	52
08	A8	Liver	Intrahepatic cholangiocarcinoma	II	F	40
09	A9	Liver	Margin of carcinoma	-	M	43
10	B1	Liver	Normal tissue	-	F	58
11	B2	Liver	Intrahepatic cholangiocarcinoma	II	F	58
12	B3	Liver	Margin of carcinoma	-	M	48
13	B4	Liver	Normal tissue	-	M	55
14	B5	Liver	Hepatocellular carcinoma	II	M	56
15	B6	Liver	Margin of carcinoma	-	M	56
16	B7	Liver	Normal tissue	-	M	65
17	B8	Liver	Clear cell carcinoma	-	M	65
18	B9	Liver	Margin of carcinoma	-	M	39
19	C1	Liver	Viral hepatitis	-	M	64
20	C2	Liver	Hepatocellular carcinoma	II	M	64
21	C3	Liver	Margin of carcinoma	-	M	64
22	C4	Liver	Normal tissue	-	M	56
23	C5	Liver	Hepatocellular carcinoma	II	M	57
24	C6	Liver	Margin of carcinoma	-	M	57
25	C7	Liver	Normal tissue	-	F	57
26	C8	Liver	Viral hepatitis	-	M	41
27	C9	Liver	Margin of carcinoma	-	M	41
28	D1	Liver	Viral hepatitis	-	M	52
29	D2	Liver	Hepatocellular carcinoma	II	M	33
30	D3	Liver	Margin of carcinoma	-	F	50
31	D4	Liver	Viral hepatitis	-	M	63
32	D5	Liver	Hepatocellular carcinoma	III	F	57
33	D6	Liver	Margin of carcinoma	-	M	46
34	D7	Liver	Adenocarcinoma	-	M	66
35	D8	Liver	Hepatocellular carcinoma	II	M	39



PROTEIN

BIOTECHNOLOGIES

36	D9	Liver	Margin of carcinoma	-	M	39
37	E1	Liver	Normal tissue	-	F	14
38	E2	Liver	Hepatocellular carcinoma	II	F	98
39	E3	Liver	Margin of carcinoma	-	M	65
40	E4	Liver	Normal tissue	-	M	56
41	E5	Liver	Hepatocellular carcinoma	III	F	70
42	E6	Liver	Margin of carcinoma	-	F	70
43	E7	Liver	Viral hepatitis	-	M	50
44	E8	Liver	Hepatocellular carcinoma	II	M	40
45	E9	Liver	Margin of carcinoma	-	M	40
46	F1	Liver	Normal tissue	-	F	1
47	F2	Liver	Hepatocellular carcinoma	I	F	69
48	F3	Liver	Margin of carcinoma	-	F	69
49	F4	Liver	Viral hepatitis	-	F	62
50	F5	Liver	Hepatocellular carcinoma	II	M	48
51	F6	Liver	Margin of carcinoma	-	M	48
52	F7	Liver	Normal tissue	-	M	63
53	F8	Liver	Hepatocellular carcinoma	II	M	46
54	F9	Liver	Margin of carcinoma	-	M	46
55	G1	Liver	Viral hepatitis	-	M	50
56	G2	Liver	Clear cell carcinoma	-	M	59
57	G3	Liver	Margin of carcinoma	-	M	58
58	G4	Liver	Normal tissue	-	M	48
59	G5	Liver	Hepatocellular carcinoma	II	M	49
60	G6	Liver	Margin of carcinoma	-	M	49
61	G7	Liver	Normal tissue	-	M	46
62	G8	Liver	Clear cell carcinoma	-	M	65
63	G9	Liver	Margin of carcinoma	-	F	40

For Research Use Only