



PROTEIN
BIOTECHNOLOGIES

HUMAN KIDNEY CANCER TISSUE ARRAY

Catalog Number: TMA-128

Each three (3) dots from three different tissue spots represents one single specimen that was selected and pathologically confirmed (Kidney carcinoma tissue array)

Cases: 23 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** CK confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Kidney	Clear cell carcinoma	-	M	54
02	A2	Kidney	Clear cell carcinoma	-	M	54
03	A3	Kidney	Clear cell carcinoma	-	M	54
04	A4	Kidney	Granular cell carcinoma	-	M	54
05	A5	Kidney	Granular cell carcinoma	-	M	54
06	A6	Kidney	Granular cell carcinoma	-	M	54
07	A7	Kidney	Clear cell carcinoma	-	M	50
08	A8	Kidney	Clear cell carcinoma	-	M	50
09	A9	Kidney	Clear cell carcinoma	-	M	50
10	B1	Kidney	Clear cell carcinoma	-	F	71
11	B2	Kidney	Clear cell carcinoma	-	F	71
12	B3	Kidney	Clear cell carcinoma	-	F	71
13	B4	Kidney	Chromophobe carcinoma	-	F	53
14	B5	Kidney	Chromophobe carcinoma	-	F	53
15	B6	Kidney	Chromophobe carcinoma	-	F	53
16	B7	Kidney	Granular cell carcinoma	-	M	67
17	B8	Kidney	Granular cell carcinoma	-	M	67
18	B9	Kidney	Granular cell carcinoma	-	M	67
19	C1	Kidney	Clear cell carcinoma	-	F	54
20	C2	Kidney	Clear cell carcinoma	-	F	54
21	C3	Kidney	Clear cell carcinoma	-	F	54
22	C4	Kidney	Clear cell carcinoma	-	F	73
23	C5	Kidney	Clear cell carcinoma	-	F	73
24	C6	Kidney	Clear cell carcinoma	-	F	73
25	C7	Kidney	Clear cell carcinoma	-	M	73
26	C8	Kidney	Connective tissue	-	M	73
27	C9	Kidney	Clear cell carcinoma	-	M	73
28	D1	Kidney	Clear cell carcinoma	-	F	64
29	D2	Kidney	Clear cell carcinoma	-	F	64
30	D3	Kidney	Clear cell carcinoma	-	F	64
31	D4	Kidney	Granular cell carcinoma	-	M	54
32	D5	Kidney	Granular cell carcinoma	-	M	54
33	D6	Kidney	Granular cell carcinoma	-	M	54
34	D7	Kidney	Clear cell carcinoma	-	M	47
35	D8	Kidney	Clear cell carcinoma	-	M	47



PROTEIN

BIOTECHNOLOGIES

36	D9	Kidney	Clear cell carcinoma	-	M	47
37	E1	Kidney	Chromophobe carcinoma	-	M	63
38	E2	Kidney	Chromophobe carcinoma	-	M	63
39	E3	Kidney	Chromophobe carcinoma	-	M	63
40	E4	Kidney	Clear cell carcinoma	-	F	50
41	E5	Kidney	Clear cell carcinoma	-	F	50
42	E6	Kidney	Clear cell carcinoma	-	F	50
43	E7	Kidney	Clear cell carcinoma	-	M	46
44	E8	Kidney	Clear cell carcinoma	-	M	46
45	E9	Kidney	Clear cell carcinoma	-	M	46
46	F1	Kidney	Clear cell carcinoma	-	M	56
47	F2	Kidney	Clear cell carcinoma	-	M	56
48	F3	Kidney	Clear cell carcinoma	-	M	56
49	F4	Kidney	Clear cell carcinoma	-	M	72
50	F5	Kidney	Clear cell carcinoma	-	M	72
51	F6	Kidney	Clear cell carcinoma	-	M	72
52	F7	Kidney	Granular cell carcinoma	-	M	51
53	F8	Kidney	Granular cell carcinoma	-	M	51
54	F9	Kidney	Granular cell carcinoma	-	M	51
55	G1	Kidney	Clear cell carcinoma	-	M	74
56	G2	Kidney	Clear cell carcinoma	-	M	74
57	G3	Kidney	Clear cell carcinoma	-	M	74
58	G4	Kidney	Clear cell carcinoma	-	M	50
59	G5	Kidney	Clear cell carcinoma	-	M	50
60	G6	Kidney	Clear cell carcinoma	-	M	50
61	G7	Kidney	Normal kidney tissue	-	F	47
62	G8	Kidney	Normal kidney tissue	-	M	23
63	G9	Kidney	Normal kidney tissue	-	M	41

For Research Use Only