



PROTEIN
BIOTECHNOLOGIES

HUMAN KIDNEY CANCER TISSUE ARRAY

Catalog Number: TMA-129

Each dot represents a DISEASED tissue spot from one individual specimen that was selected and pathologically confirmed (Kidney carcinoma tissue array)

Cases: 71 **Cores:** 72 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** CK confirmed

Layout: 8 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Kidney	Clear cell carcinoma	-	M	72
02	A2	Kidney	Granular cell carcinoma	-	M	54
03	A3	Kidney	Granular cell carcinoma	-	F	75
04	A4	Kidney	Clear cell carcinoma	-	M	62
05	A5	Kidney	Clear cell carcinoma	-	M	54
06	A6	Kidney	Transitional cell carcinoma	III	M	44
07	A7	Kidney	Clear cell carcinoma	-	F	50
08	A8	Kidney	Clear cell carcinoma	-	M	73
09	A9	Kidney	Clear cell carcinoma	-	M	50
10	B1	Kidney	Clear cell carcinoma	-	M	67
11	B2	Kidney	Transitional cell carcinoma	II	M	76
12	B3	Kidney	Clear cell carcinoma	-	M	56
13	B4	Kidney	Clear cell carcinoma	-	F	53
14	B5	Kidney	Clear cell carcinoma	-	M	45
15	B6	Kidney	Clear cell carcinoma	-	M	51
16	B7	Kidney	Granular cell carcinoma	-	M	63
17	B8	Kidney	Clear cell carcinoma	-	F	64
18	B9	Kidney	Clear cell carcinoma	-	F	71
19	C1	Kidney	Clear cell carcinoma	-	M	66
20	C2	Kidney	Granular cell carcinoma	-	F	57
21	C3	Kidney	Clear cell carcinoma	-	F	54
22	C4	Kidney	Clear cell carcinoma	-	M	63
23	C5	Kidney	Clear cell carcinoma	-	M	51
24	C6	Kidney	Transitional cell carcinoma	III	M	33
25	C7	Kidney	Clear cell carcinoma	-	M	74
26	C8	Kidney	Granular cell carcinoma	-	M	38
27	C9	Kidney	Clear cell carcinoma	-	M	54
28	D1	Kidney	Fibrous tissue and blood vessel	-	M	47
29	D2	Kidney	Clear cell carcinoma	-	F	74
30	D3	Kidney	Clear cell carcinoma	-	F	73
31	D4	Kidney	Chronic inflammation of kidney tissue	-	M	36
32	D5	Kidney	Clear cell carcinoma	-	F	64
33	D6	Kidney	Granular cell carcinoma	-	M	75
34	D7	Kidney	Clear cell carcinoma	-	M	63
35	D8	Kidney	Clear cell carcinoma	-	M	52
36	D9	Kidney	Clear cell carcinoma	-	M	46



PROTEIN

BIOTECHNOLOGIES

37	E1	Kidney	Clear cell carcinoma	-	F	50
38	E2	Kidney	Granular cell carcinoma	-	M	68
39	E3	Kidney	Clear cell carcinoma	-	M	51
40	E4	Kidney	Granular cell carcinoma	-	M	67
41	E5	Kidney	Clear cell carcinoma	-	M	60
42	E6	Kidney	Clear cell carcinoma	-	M	78
43	E7	Kidney	Granular cell carcinoma	-	M	62
44	E8	Kidney	Clear cell carcinoma	-	M	66
45	E9	Kidney	Clear cell carcinoma	-	M	65
46	F1	Kidney	Clear cell carcinoma	-	M	67
47	F2	Kidney	Clear cell carcinoma	-	F	68
48	F3	Kidney	Clear cell carcinoma	-	M	59
49	F4	Kidney	Granular cell carcinoma	-	M	61
50	F5	Kidney	Clear cell carcinoma	-	F	61
51	F6	Kidney	Clear cell carcinoma	-	M	62
52	F7	Kidney	Clear cell carcinoma	-	M	50
53	F8	Kidney	Clear cell carcinoma	-	M	51
54	F9	Kidney	Clear cell carcinoma	-	M	70
55	G1	Kidney	Transitional cell carcinoma	III	M	71
56	G2	Kidney	Clear cell carcinoma	-	F	61
57	G3	Kidney	Granular cell carcinoma	-	F	47
58	G4	Kidney	Clear cell carcinoma	-	M	52
59	G5	Kidney	Clear cell carcinoma	-	M	72
60	G6	Kidney	Clear cell carcinoma	-	M	58
61	G7	Kidney	Clear cell carcinoma	-	F	58
62	G8	Kidney	Granular cell carcinoma	-	F	33
63	G9	Kidney	Granular cell carcinoma	-	F	45
64	H1	Kidney	Clear cell carcinoma	-	M	53
65	H2	Kidney	Transitional cell carcinoma	II	M	68
66	H3	Kidney	Clear cell carcinoma	-	M	45
67	H4	Kidney	Transitional cell carcinoma	III	M	62
68	H5	Kidney	Clear cell carcinoma	-	M	59
69	H6	Kidney	Clear cell carcinoma	-	M	59
70	H7	Kidney	Normal kidney tissue	-	M	20
71	H8	Kidney	Normal kidney tissue	-	M	71
72	H9	Kidney	Normal kidney tissue	-	M	70

For Research Use Only