



PROTEIN
BIOTECHNOLOGIES

HUMAN KIDNEY CANCER TISSUE ARRAY

Catalog Number: TMA-130

Each three (3) dots from three different tissue spots represents one single specimen that was selected and pathologically confirmed (Normal kidney tissue array)

Cases: 21 **Cores:** 54 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** CK confirmed

Layout: 6 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Kidney	Cortex	-	M	51
02	A2	Kidney	Cortex	-	M	51
03	A3	Kidney	Cortex	-	M	51
04	A4	Kidney	Cortex	-	M	41
05	A5	Kidney	Cortex	-	M	41
06	A6	Kidney	Cortex	-	M	41
07	A7	Kidney	Cortex	-	M	59
08	A8	Kidney	Cortex	-	M	59
09	A9	Kidney	Cortex	-	M	59
10	B1	Kidney	Cortex	-	M	65
11	B2	Kidney	Cortex	-	M	65
12	B3	Kidney	Cortex	-	M	65
13	B4	Kidney	Cortex	-	M	72
14	B5	Kidney	Cortex	-	M	72
15	B6	Kidney	Cortex	-	M	72
16	B7	Kidney	Cortex	-	M	52
17	B8	Kidney	Cortex	-	M	52
18	B9	Kidney	Cortex	-	M	52
19	C1	Kidney	Cortex	-	F	58
20	C2	Kidney	Cortex	-	F	58
21	C3	Kidney	Cortex	-	F	58
22	C4	Kidney	Cortex	-	F	45
23	C5	Kidney	Cortex	-	F	45
24	C6	Kidney	Cortex	-	F	45
25	C7	Kidney	Cortex	-	M	60
26	C8	Kidney	Cortex	-	M	60
27	C9	Kidney	Cortex	-	M	60
28	D1	Kidney	Cortex	-	M	41
29	D2	Kidney	Cortex	-	M	41
30	D3	Kidney	Cortex	-	M	41
31	D4	Kidney	Cortex	-	M	63
32	D5	Kidney	Cortex	-	M	63
33	D6	Kidney	Cortex	-	M	63
34	D7	Kidney	Medulla	-	M	23
35	D8	Kidney	Cortex	-	M	23
36	D9	Kidney	Medulla	-	M	23



PROTEIN

BIOTECHNOLOGIES

37	E1	Kidney	Cortex	-	F	51
38	E2	Kidney	Cortex	-	F	51
39	E3	Kidney	Cortex	-	F	51
40	E4	Kidney	Cortex	-	M	18
41	E5	Kidney	Cortex	-	M	18
42	E6	Kidney	Cortex	-	M	18
43	E7	Kidney	Cortex	-	F	47
44	E8	Kidney	Cortex	-	F	47
45	E9	Kidney	Cortex	-	F	47
46	F1	Kidney	Cortex	-	M	76
47	F2	Kidney	Cortex	-	M	76
48	F3	Kidney	Cortex	-	M	76
49	F4	Kidney	Cortex	-	F	53
50	F5	Kidney	Cortex	-	F	53
51	F6	Kidney	Cortex	-	F	53
52	F7	Kidney	Clear cell carcinoma	-	M	54
53	F8	Kidney	Clear cell carcinoma	-	M	63
54	F9	Kidney	Granulosa clear cell carcinoma	-	M	45

For Research Use Only