



**PROTEIN**  
BIOTECHNOLOGIES

**HUMAN THYROID CANCER TISSUE ARRAY**

**Catalog Number:** TMA-136

Each three (3) dots from three different tissue spots represents one single specimen that was selected and pathologically confirmed

**Cases:** 63    **Cores:** 63    **Diameter:** 1.5mm    **Thickness:** 5 µm    **Standard IHC:** Tg confirmed

**Layout:** 7 x 9

	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>
<b>A</b>	1	2	3	4	5	6	7	8	9
<b>B</b>	10	11	12	13	14	15	16	17	18
<b>C</b>	19	20	21	22	23	24	25	26	27
<b>D</b>	28	29	30	31	32	33	34	35	36
<b>E</b>	37	38	39	40	41	42	43	44	45
<b>F</b>	46	47	48	49	50	51	52	53	54
<b>G</b>	55	56	57	58	59	60	61	62	63

<b>Core No.</b>	<b>Location</b>	<b>Tissue</b>	<b>Diagnosis</b>	<b>Grade</b>	<b>Sex</b>	<b>Age</b>
01	A1	Thyroid	Follicular carcinoma	-	F	70
02	A2	Thyroid	Follicular carcinoma	-	F	70
03	A3	Thyroid	Follicular carcinoma	-	F	70
04	A4	Thyroid	Papillary carcinoma	-	M	46
05	A5	Thyroid	Papillary carcinoma	-	M	46
06	A6	Thyroid	Papillary carcinoma	-	M	46
07	A7	Thyroid	Squamous cell of carcinoma	I	M	45
08	A8	Thyroid	Squamous cell of carcinoma	I	M	45
09	A9	Thyroid	Squamous cell of carcinoma	I	M	45
10	B1	Thyroid	Follicular variant papillary carcinoma	-	M	50
11	B2	Thyroid	Follicular variant papillary carcinoma	-	M	50
12	B3	Thyroid	Follicular variant papillary carcinoma	-	M	50
13	B4	Thyroid	Papillary carcinoma	-	F	68
14	B5	Thyroid	Papillary carcinoma	-	F	68
15	B6	Thyroid	Papillary carcinoma	-	F	68
16	B7	Thyroid	Papillary carcinoma	-	F	28
17	B8	Thyroid	Papillary carcinoma	-	F	28
18	B9	Thyroid	Papillary carcinoma	-	F	28
19	C1	Thyroid	Papillary carcinoma	-	M	27
20	C2	Thyroid	Papillary carcinoma	-	M	27
21	C3	Thyroid	Papillary carcinoma	-	M	27
22	C4	Thyroid	Follicular variant papillary carcinoma	-	F	49
23	C5	Thyroid	Follicular variant papillary carcinoma	-	F	49
24	C6	Thyroid	Follicular variant papillary carcinoma	-	F	49
25	C7	Thyroid	Follicular variant papillary carcinoma	-	F	50
26	C8	Thyroid	Follicular variant papillary carcinoma	-	F	50
27	C9	Thyroid	Fibrous tissue and striated muscles	-	F	50
28	D1	Thyroid	Papillary carcinoma	-	F	48
29	D2	Thyroid	Papillary carcinoma	-	F	48
30	D3	Thyroid	Papillary carcinoma	-	F	48
31	D4	Thyroid	Papillary carcinoma	-	F	45
32	D5	Thyroid	Papillary carcinoma	-	F	45
33	D6	Thyroid	Papillary carcinoma	-	F	45
34	D7	Thyroid	Papillary carcinoma	-	M	40
35	D8	Thyroid	Papillary carcinoma	-	M	40



# PROTEIN

BIOTECHNOLOGIES

36	D9	Thyroid	Papillary carcinoma	-	M	40
37	E1	Thyroid	Follicular variant papillary carcinoma	-	F	25
38	E2	Thyroid	Follicular variant papillary carcinoma	-	F	25
39	E3	Thyroid	Follicular variant papillary carcinoma	-	F	25
40	E4	Thyroid	Papillary carcinoma	-	F	45
41	E5	Thyroid	Papillary carcinoma	-	F	45
42	E6	Thyroid	Papillary carcinoma	-	F	45
43	E7	Thyroid	Papillary carcinoma	-	F	41
44	E8	Thyroid	Papillary carcinoma	-	F	41
45	E9	Thyroid	Papillary carcinoma	-	F	41
46	F1	Thyroid	Papillary carcinoma	-	F	46
47	F2	Thyroid	Papillary carcinoma	-	F	46
48	F3	Thyroid	Papillary carcinoma	-	F	46
49	F4	Thyroid	Papillary carcinoma	-	F	16
50	F5	Thyroid	Papillary carcinoma	-	F	16
51	F6	Thyroid	Papillary carcinoma	-	F	16
52	F7	Thyroid	Undifferentiated carcinoma	-	F	81
53	F8	Thyroid	Undifferentiated carcinoma	-	F	81
54	F9	Thyroid	Undifferentiated carcinoma	-	F	81
55	G1	Thyroid	Papillary carcinoma	-	F	34
56	G2	Thyroid	Papillary carcinoma	-	F	34
57	G3	Thyroid	Papillary carcinoma	-	F	34
58	G4	Thyroid	Papillary carcinoma	-	F	60
59	G5	Thyroid	Papillary carcinoma	-	F	60
60	G6	Thyroid	Papillary carcinoma	-	F	60
61	G7	Thyroid	Normal thyroid tissue	-	F	27
62	G8	Thyroid	Normal thyroid tissue	-	F	40
63	G9	Thyroid	Normal thyroid tissue	-	F	64

**For Research Use Only**