



PROTEIN BIOTECHNOLOGIES

HUMAN LYMPHOMA CANCER TISSUE ARRAY

Catalog Number: TMA-147

Each dot represents a DISEASED tissue spot from one individual specimen that was selected and pathologically confirmed

Cases: 60 **Cores:** 60 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** LCA confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Lymph node	Reactive hyperplasia of lymph nodes	-	F	46
02	A2	Lymph node	Plasmacytoid lymphocytic lymphoma	-	M	50
03	A3	Intestine	B cell mucosa-associated lymphoma	-	M	35
04	A4	Lymph node	Large B cell lymphoma	-	M	42
05	A5	Lymph node	Large B cell lymphoma of brain	-	F	46
06	A6	Lymph node	Small B cell lymphoma	-	F	36
07	A7	Lymph node	B cell lymphoma	-	M	53
08	A8	Lymph node	Large B cell lymphoma	-	M	54
09	A9	Lymph node	B cell lymphoma	-	F	44
10	B1	Lymph node	Burkitt's lymphoma	-	M	7
11	B2	Lymph node	Hodgkin lymphoma (nodular sclerosis type)	-	M	17
12	B3	Lymph node	Angioimmunoblastic lymphadenopathic T cell lymphoma	-	M	67
13	B4	Lymph node	T lymphoblastic lymphoma	-	M	8
14	B5	Lymph node	Chronic inflammation of lymph node	-	M	50
15	B6	Lymph node	Diffuse small cleaved lymphoma	-	M	53
16	B7	Lymph node	Large B cell lymphoma	-	M	51
17	B8	Lymph node	Plasmacytic lymphoma of nuchae	-	F	60
18	B9	Lymph node	B cell lymphoma	-	M	42
19	C1	Lymph node	Malignant lymphoma (No T cell and No B cell)	-	M	64
20	C2	Lymph node	B cell lymphoma	-	M	64
21	C3	Lymph node	B cell lymphoma	-	M	64
22	C4	Lymph node	B cell lymphoma	-	F	50
23	C5	Lymph node	Diffuse large B cell lymphoma	-	M	55
24	C6	Kidney	Plasmacytoid T cell lymphoma	-	F	58
25	C7	Lymph node	B cell mucosa-associated lymphoma of small intestine	-	M	41
26	C8	Lymph node	B cell lymphoma	-	M	35
27	C9	Lymph node	B cell lymphoma	-	M	71
28	D1	Lymph node	Hodgkin lymphoma (mixed type)	-	M	4
29	D2	Lymph node	B cell lymphoma	-	F	71
30	D3	Lymph node	B cell lymphoma	-	M	45
31	D4	Lymph node	B cell lymphoma	-	F	44
32	D5	Lymph node	B cell lymphoma	-	M	67
33	D6	Lymph node	Large B cell lymphoma	-	F	46



PROTEIN BIOTECHNOLOGIES

34	D7	Lymph node	B cell lymphoma	-	M	58
35	D8	Lymph node	B cell mucosa-associated lymphoma	-	M	81
36	D9	Lymph node	Diffuse B cell lymphoma	-	M	27
37	E1	Lymph node	Diffuse B cell lymphoma	-	F	63
38	E2	Lymph node	B cell lymphoma	-	M	51
39	E3	Lymph node	B cell lymphoma	-	M	15
40	E4	Tonsil	B cell lymphoma	-	M	38
41	E5	Colon	B cell lymphoma of right hemicolon	-	M	50
42	E6	Lung	Large B cell lymphoma	-	F	68
43	E7	Intestine	B cell mucosa-associated lymphoma	-	F	36
44	E8	Stomach	B cell mucosa-associated lymphoma of stomach	-	M	43
45	E9	Lymph node	B cell lymphoma	-	M	21
46	F1	Lymph node	Diffuse large B cell lymphoma	-	M	68
47	F2	Lymph node	B cell lymphoma	-	F	63
48	F3	Spleen	Large B cell lymphoma	-	M	27
49	F4	Lymph node	Anaplastic large cell lymphoma	-	M	69
50	F5	Lymph node	B cell lymphoma	-	M	65
51	F6	Lymph node	B cell lymphoma	-	M	45
52	F7	Lymph node	Centroblastic lymphoma	-	F	76
53	F8	Stomach	B cell lymphoma	-	M	76
54	F9	Kidney	Small B cell Lymphoma	-	F	54
55	G1	Lymph node	B cell lymphoma	-	F	66
56	G2	Lymph node	Diffuse B cell lymphoma	-	M	13
57	G3	Lymph node	Diffuse large cell lymphoma	-	F	64
58	G4	Lymph node	Hodgkin lymphoma (lymphocyte predominance type)	-	F	73
59	G5	Lymph node	Mantle cell lymphoma	-	M	57
60	G6	Lymph node	Follicular center cell lymphoma		F	68

For Research Use Only