



PROTEIN
BIOTECHNOLOGIES

HUMAN HEAD AND NECK CANCER TISSUE ARRAY

Catalog Number: TMA-152

Each three (3) dots from three different tissue spots represents one single specimen that was selected and pathologically confirmed

Cases: 24 **Cores:** 63 **Diameter:** 1.5mm **Thickness:** 5 µm **Standard IHC:** H-CK confirmed

Layout: 7 x 9

	1	2	3	4	5	6	7	8	9
A	1	2	3	4	5	6	7	8	9
B	10	11	12	13	14	15	16	17	18
C	19	20	21	22	23	24	25	26	27
D	28	29	30	31	32	33	34	35	36
E	37	38	39	40	41	42	43	44	45
F	46	47	48	49	50	51	52	53	54
G	55	56	57	58	59	60	61	62	63

Core No.	Location	Tissue	Diagnosis	Grade	Sex	Age
01	A1	Cheek part	Squamous cell carcinoma of cheek part (grade I)	I	F	51
02	A2	Cheek part	Squamous cell carcinoma of cheek part (grade II)	II	F	51
03	A3	Cheek part	Squamous cell carcinoma of cheek part (grade I)	I	F	51
04	A4	Tongue	Squamous cell carcinoma of tongue (grade II)	II	M	53
05	A5	Tongue	Squamous cell carcinoma of tongue (grade II)	II	M	53
06	A6	Tongue	Squamous cell carcinoma of tongue (grade II)	II	M	53
07	A7	Larynx	Squamous cell carcinoma of larynx (grade II)	II	M	65
08	A8	Larynx	Squamous cell carcinoma of larynx (grade II)	II	M	65
09	A9	Larynx	Squamous cell carcinoma of larynx (grade II)	II	M	65
10	B1	Cheek part	Squamous cell carcinoma of cheek part (grade I)	I	M	52
11	B2	Cheek part	Squamous cell carcinoma of cheek part (grade I)	I	M	52
12	B3	Cheek part	Squamous cell carcinoma of cheek part (grade II)	II	M	52
13	B4	Tongue	Squamous cell carcinoma of tongue (grade III)	III	M	60
14	B5	Tongue	Squamous cell carcinoma of tongue (grade III)	III	M	60
15	B6	Tongue	Squamous cell carcinoma of tongue (grade III)	III	M	60
16	B7	Larynx	Squamous cell carcinoma of larynx (grade III)	III	M	52
17	B8	Larynx	Squamous cell carcinoma of larynx (grade III)	III	M	52
18	B9	Larynx	Squamous cell carcinoma of larynx (grade III)	III	M	52
19	C1	Cheek part	Squamous cell carcinoma of cheek part (grade I)	I	M	49
20	C2	Cheek part	Squamous cell carcinoma of cheek part (grade I)	I	M	49
21	C3	Cheek part	Squamous cell carcinoma of cheek part (grade I)	I	M	49
22	C4	Tongue	Squamous cell carcinoma of tongue (grade II)	II	M	57
23	C5	Tongue	Squamous cell carcinoma of tongue (grade II)	II	M	57
24	C6	Tongue	Squamous cell carcinoma of tongue (grade II)	II	M	57
25	C7	Larynx	Squamous cell carcinoma of larynx (grade I)	I	M	56
26	C8	Larynx	Squamous cell carcinoma of larynx (grade I)	I	M	56
27	C9	Larynx	Squamous cell carcinoma of larynx (grade I)	I	M	56
28	D1	Cheek part	Squamous cell carcinoma of cheek part (grade III)	III	F	52
29	D2	Cheek part	Squamous cell carcinoma of cheek part (grade III)	III	F	52
30	D3	Cheek part	Squamous cell carcinoma of cheek part (grade III)	III	F	52
31	D4	Tongue	Squamous cell carcinoma of tongue (grade I)	I	F	48
32	D5	Tongue	Squamous cell carcinoma of tongue (grade I)	I	F	48
33	D6	Tongue	Squamous cell carcinoma of tongue (grade I)	I	F	48
34	D7	Larynx	Squamous cell carcinoma of larynx (grade III)	III	M	46
35	D8	Larynx	Squamous cell carcinoma of larynx (grade III)	III	M	46



PROTEIN

BIOTECHNOLOGIES

36	D9	Larynx	Squamous cell carcinoma of larynx (grade III)	III	M	46
37	E1	Upper jaw	Squamous cell carcinoma of upper jaw (grade II)	II	M	54
38	E2	Upper jaw	Squamous cell carcinoma of upper jaw (grade II)	II	M	54
39	E3	Upper jaw	Squamous cell carcinoma of upper jaw (grade II)	II	M	54
40	E4	Tongue	Squamous cell carcinoma of tongue (grade III)	III	M	50
41	E5	Tongue	Squamous cell carcinoma of tongue (grade III)	III	M	50
42	E6	Tongue	Squamous cell carcinoma of tongue (grade II)	II	M	50
43	E7	Larynx	Squamous cell carcinoma of larynx (grade II)	II	F	69
44	E8	Larynx	Squamous cell carcinoma of larynx (grade II)	II	F	69
45	E9	Larynx	Squamous cell carcinoma of larynx (grade II)	II	F	69
46	F1	Larynx	Squamous cell carcinoma of larynx (grade II)	II	M	51
47	F2	Larynx	Squamous cell carcinoma of larynx (grade II)	II	M	51
48	F3	Larynx	Squamous cell carcinoma of larynx (grade II)	II	M	51
49	F4	Nose	Squamous cell carcinoma of nasopharynx (grade III))	M	50
50	F5	Nose	Squamous cell carcinoma of nasopharynx (grade III))	M	50
51	F6	Nose	Squamous cell carcinoma of nasopharynx (grade III))	M	50
52	F7	Nose	Squamous cell carcinoma of nasopharynx (grade III))	M	62
53	F8	Nose	Squamous cell carcinoma of nasopharynx (grade III))	M	62
54	F9	Nose	Squamous cell carcinoma of nasopharynx (grade III))	M	62
55	G1	Nose	Squamous cell carcinoma of nasopharynx (grade III))	M	50
56	G2	Nose	Squamous cell carcinoma of nasopharynx (grade III))	M	50
57	G3	Nose	Squamous cell carcinoma of nasopharynx (grade III))	M	50
58	G4	Nose	Squamous cell carcinoma of nasopharynx (grade III))	F	56
59	G5	Nose	Squamous cell carcinoma of nasopharynx (grade III))	F	56
60	G6	Nose	Squamous cell carcinoma of nasopharynx (grade III))	F	56
61	G7	Gingiva	Mucitis of gingiva	-	M	33
62	G8	Gingiva	Mucosa tissues of gingiva	-	M	45
63	G9	Gingiva	Mucitis of gingiva	-	M	41

For Research Use Only