



**PROTEIN**  
BIOTECHNOLOGIES

**HUMAN HEAD AND NECK CANCER TISSUE ARRAY**

**Catalog Number:** TMA-153

Each three (3) cores from three different tissue spots represents one single specimen that was pathologically selected and confirmed

**Cases:** 72    **Cores:** 72    **Diameter:** 1.5mm    **Thickness:** 5 µm    **Standard IHC:** H-CK confirmed

**Layout:** 8 x 9

	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>
<b>A</b>	1	2	3	4	5	6	7	8	9
<b>B</b>	10	11	12	13	14	15	16	17	18
<b>C</b>	19	20	21	22	23	24	25	26	27
<b>D</b>	28	29	30	31	32	33	34	35	36
<b>E</b>	37	38	39	40	41	42	43	44	45
<b>F</b>	46	47	48	49	50	51	52	53	54
<b>G</b>	55	56	57	58	59	60	61	62	63
<b>H</b>	64	65	66	67	68	69	70	71	72

<b>Core No.</b>	<b>Location</b>	<b>Tissue</b>	<b>Diagnosis</b>	<b>Grade</b>	<b>Sex</b>	<b>Age</b>
01	A1	Cheek part	Squamous cell carcinoma	II	F	51
02	A2	Cheek part	Squamous cell carcinoma	II	F	51
03	A3	Cheek part	Squamous cell carcinoma	II	F	51
04	A4	Upper jaw	Squamous cell carcinoma	I	M	40
05	A5	Upper jaw	Squamous cell carcinoma	I	M	40
06	A6	Upper jaw	Squamous cell carcinoma	II	M	40
07	A7	Larynx	Squamous cell carcinoma	III	M	65
08	A8	Larynx	Squamous cell carcinoma	III	M	65
09	A9	Larynx	Squamous cell carcinoma	III	M	65
10	B1	Cheek part	Squamous cell carcinoma	I	M	52
11	B2	Cheek part	Squamous cell carcinoma	I	M	52
12	B3	Cheek part	Squamous cell carcinoma	I	M	52
13	B4	Tongue	Squamous cell carcinoma	III	M	60
14	B5	Tongue	Squamous cell carcinoma	III	M	60
15	B6	Tongue	Squamous cell carcinoma	III	M	60
16	B7	Larynx	Squamous cell carcinoma	III	M	52
17	B8	Larynx	Squamous cell carcinoma	III	M	52
18	B9	Larynx	Squamous cell carcinoma	III	M	52
19	C1	lip	Squamous cell carcinoma	III	M	57
20	C2	lip	Squamous cell carcinoma	III	M	57
21	C3	lip	Squamous cell carcinoma	III	M	57
22	C4	Tongue	Squamous cell carcinoma	II	M	57
23	C5	Tongue	Squamous cell carcinoma	II	M	57
24	C6	Tongue	Squamous cell carcinoma	II	M	57
25	C7	Larynx	Squamous cell carcinoma	I	M	56
26	C8	Larynx	Squamous cell carcinoma	I	M	56
27	C9	Larynx	Squamous cell carcinoma	I	M	56
28	D1	lip	Squamous cell carcinoma	II	M	81
29	D2	lip	Squamous cell carcinoma	II	M	81
30	D3	lip	Squamous cell carcinoma	II	M	81
31	D4	Tongue	Squamous cell carcinoma	III	F	48
32	D5	Tongue	Squamous cell carcinoma	III	F	48
33	D6	Tongue	Squamous cell carcinoma	III	F	48



# PROTEIN

BIOTECHNOLOGIES

34	D7	Larynx	Squamous cell carcinoma	III	M	46
35	D8	Larynx	Squamous cell carcinoma	III	M	46
36	D9	Larynx	Squamous cell carcinoma	III	M	46
37	E1	Gingiva	Squamous cell carcinoma	II	M	49
38	E2	Gingiva	Squamous cell carcinoma	II	M	49
39	E3	Gingiva	Squamous cell carcinoma	II	M	49
40	E4	Tongue	Squamous cell carcinoma	III	M	50
41	E5	Tongue	Squamous cell carcinoma	III	M	50
42	E6	Tongue	Squamous cell carcinoma	III	M	50
43	E7	Larynx	Squamous cell carcinoma	II	F	69
44	E8	Larynx	Squamous cell carcinoma	II	F	69
45	E9	Larynx	Squamous cell carcinoma	II	F	69
46	F1	Nose	Squamous cell carcinoma	III	F	63
47	F2	Nose	Squamous cell carcinoma	III	F	63
48	F3	Nose	Squamous cell carcinoma	III	F	63
49	F4	Nose	Squamous cell carcinoma	III	M	43
50	F5	Nose	Squamous cell carcinoma	III	M	43
51	F6	Nose	Squamous cell carcinoma	III	M	43
52	F7	Nasopharynx	Squamous cell carcinoma	III	M	62
53	F8	Nasopharynx	Squamous cell carcinoma	III	M	62
54	F9	Nasopharynx	Squamous cell carcinoma	III	M	62
55	G1	Nose	Squamous cell carcinoma	II	M	42
56	G2	Nose	Squamous cell carcinoma	II	M	42
57	G3	Nose	Squamous cell carcinoma	II	M	42
58	G4	Ethmoid sinus	Squamous cell carcinoma	III	M	34
59	G5	Ethmoid sinus	Squamous cell carcinoma	III	M	34
60	G6	Ethmoid sinus	Squamous cell carcinoma	III	M	34
61	G7	Gingiva	Mucosa	-	M	33
62	G8	Gingiva	Mucosa	-	M	33
63	G9	Gingiva	Mucosa	-	M	45
64	H1	Gingiva	Mucosa	-	M	45
65	H2	Gingiva	Mucosa	-	M	41
66	H3	Gingiva	Mucosa	-	M	41
67	H4	Tongue	Mucosa	-	F	51
68	H5	Tongue	Mucosa	-	F	51
69	H6	Tongue	Mucosa	-	F	39
70	H7	Tongue	Mucosa	-	F	39
71	H8	Tongue	Mucosa	-	M	66
72	H9	Tongue	Mucosa	-	M	66